

Long-read pharmacogenomics (PGx) with StarPhase

StarPhase is a PGx diplotyping tool that leverages long-read sequencing to resolve complex pharmacogenes. Key features include:

- **Diploypes 21 CPIC Class A pharmacogenes** – Including complex genes *HLA-A*, *HLA-B*, and *CYP2D6*
- **Phase-aware diplotyping** – Generates haplotype-resolved PGx diploypes using long reads
- **Full-length consensus sequences** – Generates complete genomic sequences for *HLA* genes and *CYP2D6*, enabling novel allele discovery at basepair resolution
- **Custom visualizations** – Creates sample-specific IGV visualizations to enable exploration of complex loci

StarPhase generates high-accuracy PGx diplotype calls




Figure 1. Comparison to HPRC benchmark. The vast majority of diploypes calls match the benchmark (98.1%). An additional 1.5% are caused by minor discrepancies or database updates. Manual inspection of the final nine diploypes (0.3%) supports the StarPhase diploypes calls.

StarPhase updates PGx diplotype calls for GeT-RM samples




Figure 2. Comparison to GeT-RM sample benchmark. The majority of diploypes calls match the GeT-RM call set (73.8%). Of the differences, 38 diploypes (14.8%) match new haplotypes that were not defined when the benchmark was created. For *UGT1A1*, 22 diploypes (8.6%) contain a retired haplotype, and are updated accordingly. Of the final seven mismatches, six can be explained by limitations of the benchmark assays. In total, StarPhase provided updated results for 26.2% of diploypes in this GeT-RM sample benchmark.

Extending the methods to nine additional HLA genes




Figure 3. HLA-DRB1 in CEPH pedigree 1463. StarPhase has been extended to report nine more complex HLA genes. This figure shows the diploypes calls for *HLA-DRB1* in CEPH pedigree 1463. StarPhase identifies each diplotype independently and generates no inheritance violations for the pedigree. Full list of additional HLA genes: *HLA-C*, *DPA1*, *DPB1*, *DQA1*, *DQB1*, *DRB1*, *DRB3*, *DRB4*, and *DRB5*.

StarPhase visualizations enhance interpretability of long-read PGx data





Figure 4. Example *4.004x2 haplotype. This IGV screenshot has been enhanced to highlight key evidence from a StarPhase *4.004x2 haplotype call. Assembled regions are labeled on the bottom track, with two *4.004 alleles highlighted in pink. Long reads that span both copies of *4.004 are highlighted in orange. These reads provide direct observational evidence of the duplication.

Conclusion

StarPhase is a long-read pharmacogenomic diplotyper that:

- Provides highly accurate diplotype results from long-read observations
- Provides refined PGx diploypes for GeT-RM benchmark samples
- Generates full-length haplotype sequences and visualizations for complex pharmacogenes



Available on GitHub:

<https://github.com/PacificBiosciences/pb-StarPhase>

References

1. Behera, Sairam, et al. "Comprehensive and accurate genome analysis at scale using DRAGEN accelerated algorithms." *bioRxiv* (2024).
2. Chen, Xiao, et al. "Cyrius: accurate CYP2D6 genotyping using whole-genome sequencing data." *The pharmacogenomics journal* 21.2 (2021): 251-261.
3. Lai, Sheng-Kai, et al. "A novel framework for human leukocyte antigen (HLA) genotyping using probe capture-based targeted next-generation sequencing and computational analysis." *Computational and Structural Biotechnology Journal* 23 (2024): 1562-1571.
4. Sangkuhl, Katrin, et al. "Pharmacogenomics clinical annotation tool (PharmCAT)." *Clinical Pharmacology & Therapeutics* 107.1 (2020): 203-210.
5. Wang, Ting, et al. "The Human Pangenome Project: a global resource to map genomic diversity." *Nature* 604.7906 (2022): 437-446.