



PACIFIC
BIOSCIENCES®



PacBio HiFiViral Workflow Overview: Multiplexed Amplicon Library Preparation for Full-Viral Genome Sequencing of SARS-CoV-2

Sequel II System ICS v9.0 / Sequel II Chemistry 2.0 / SMRT Link v9.0

Sequel IIe System ICS v10.0 / Sequel II Chemistry 2.0 / SMRT Link v10.0

For Research Use Only. Not for use in diagnostic procedures. © Copyright 2021 by Pacific Biosciences of California, Inc. All rights reserved.

PN 102-084-800 Version 2021-02-25-A (February 2021)

Multiplexed Amplicon Library Preparation for Full-Viral Genome Sequencing of SARS-CoV-2

1. PacBio HiFiViral for SARS-CoV-2 Workflow Overview
2. Multiplexed SARS-CoV-2 Amplicon Library Preparation Using PacBio Barcoded M13 Primers
3. Multiplexed SARS-CoV-2 Amplicon Library Sequencing Workflow Recommendations
4. Multiplexed SARS-CoV-2 Amplicon Data Analysis Recommendations
5. Technical Documentation & Applications Support Resources

Appendix: RNA Isolation Kit Options for Full-Viral Genome Sequencing of SARS-CoV-2

SARS-CoV-2 FULL-VIRAL GENOME SEQUENCING: HOW TO GET STARTED

Application-Specific Resources Website

Application-Specific Procedure & Checklist

Application Consumable Bundle Purchasing Guide

Library Construction, Sequencing & Analysis

COVID-19 SEQUENCING TOOLS AND RESOURCES

UNDERSTAND THE EVOLVING PANDEMIC WITH PACBIO SEQUENCING

Targeted Sequencing for SARS-CoV-2 Excellence

1.2 kb amplicons

20 kb

Figure 1: 20 x 1.2 kb amplicons tiled across the SARS-CoV-2 genome.

This procedure involves two rounds of amplification, first using M13-tailed target-specific primers to tail the PCR products with a universal M13 sequence, followed by a second round PCR using barcoded M13 primers. There are 32 forward (F) and 32 reverse (R) barcoded M13 primers available for barcoding and when used in different combinations allow multiplexing up to 1024 samples.

The number of samples that can be multiplexed per SMRT Cell depends on the desired coverage per sample. For Sequel II and IIe Systems, PacBio recommends multiplexing up to 900-900 SARS-CoV-2 samples on a single SMRT Cell SM.

For any questions or additional information about this procedure, please contact support@pacb.com.

Page 1 Part Number 102-075-000 Version 03 (February 11, 2021)
SARS-CoV-2 Full-Genome Sequencing Workflow - Collaboration between LACOG and Pacific

[PacBio COVID-19 Sequencing Tools and Resources Website](#)

Summary overview of application-specific sample preparation and data analysis workflow recommendations

Customer Collaboration – PacBio HiFiViral for SARS-CoV-2 Workflow: High-Throughput Multiplexing 1.2 kb Amplicons for Full-Viral Genome Sequencing of SARS-CoV-2

This protocol is the result of a collaboration between LACOG and Pacific Biosciences. We acknowledge the contributions of:

- LACOG: Brian Kivinger, Michael Levandouki, Jonathan Williams, Brian Novell, Scott Palmer, Ashley Letovsky, Gian Zeng, Xiaoyi Maria Estenberg
- Pacific Biosciences: Peter Bockarie, George Yuan, Elizabeth Thong, Jonas Kotlach

This document describes a workflow for whole viral genome sequencing of SARS-CoV-2 samples on the Sequel II and Sequel IIe Systems using a targeted PCR approach. PCR strand cDNA is synthesized using the SuperScript IV (S-IV) cDNA synthesis kit (Thermo Fisher Scientific). The resulting cDNA is then PCR amplified, asymmetrically barcoded and subsequently pooled with other samples for SMRTbell library construction and multiplex sequencing on the Sequel II and IIe Systems.

This procedure recommends amplification and sequencing of twenty nine overlapping 1.2 kb amplicons that are tiled across the full 29.8 kb SARS-CoV-2 genome (Figure 1). The PCR primers were designed by Blake E. Freed of the School of Natural and Computational Sciences, Massey University, Auckland, New Zealand.

1.2 kb amplicons

20 kb

Figure 1: 20 x 1.2 kb amplicons tiled across the SARS-CoV-2 genome.

This procedure involves two rounds of amplification, first using M13-tailed target-specific primers to tail the PCR products with a universal M13 sequence, followed by a second round PCR using barcoded M13 primers. There are 32 forward (F) and 32 reverse (R) barcoded M13 primers available for barcoding and when used in different combinations allow multiplexing up to 1024 samples.

The number of samples that can be multiplexed per SMRT Cell depends on the desired coverage per sample. For Sequel II and IIe Systems, PacBio recommends multiplexing up to 900-900 SARS-CoV-2 samples on a single SMRT Cell SM.

For any questions or additional information about this procedure, please contact support@pacb.com.

Page 1 Part Number 102-075-000 Version 03 (February 11, 2021)
SARS-CoV-2 Full-Genome Sequencing Workflow - Collaboration between LACOG and Pacific

[Customer Collaboration – PacBio HiFiViral for SARS-CoV-2 Workflow: High-Throughput Multiplexing 1.2 kb Amplicons for Full-Viral Genome Sequencing of SARS-CoV-2](#) (102-075-000 / 102-082-500)

Technical documentation containing sample library construction and sequencing preparation protocol details

Application Consumable Bundles
Generate Highly Accurate Long-Read Sequencing Data You Can Trust

With this PacBio® Application Consumable Purchasing Guide, you can easily order the required consumables for the Sequel® II System. Simply choose your SMRT® Sequencing Application and with the single part number price you order to get started!

Application	Name and Part Number	# of Samples	Contents and Quantities
HiFi Reads for de novo Assembly and Variant Detection	Sequel II HiFi Reads-10 (101-900-000)	10	1000x PacBio Sequencing Program 2.0 (101-900-000) 10x Sequel II HiFi Read Cell (101-900-000) 10x Sequel II HiFi Read Cell (101-900-000) 10x Sequel II HiFi Read Cell (101-900-000)
De novo Assembly for Low DNA Input Samples	Sequel II De Novo Low DNA Input Kit (101-900-000)	10	1000x PacBio Sequencing Program 2.0 (101-900-000) 10x Sequel II HiFi Read Cell (101-900-000) 10x Sequel II HiFi Read Cell (101-900-000) 10x Sequel II HiFi Read Cell (101-900-000)
De novo Assembly for Research Applications	Sequel II Multiplexed Assembly Kit-10 (101-900-000)	10	1000x PacBio Sequencing Program 2.0 (101-900-000) 10x Sequel II HiFi Read Cell (101-900-000) 10x Sequel II HiFi Read Cell (101-900-000) 10x Sequel II HiFi Read Cell (101-900-000)
Structural Variant Detection	Sequel II Multiplexed SV Detection Kit-10 (101-900-000)	10	1000x PacBio Sequencing Program 2.0 (101-900-000) 10x Sequel II HiFi Read Cell (101-900-000) 10x Sequel II HiFi Read Cell (101-900-000) 10x Sequel II HiFi Read Cell (101-900-000)
De novo Assembly for Phasing	Sequel II De Novo Phasing Kit-10 (101-900-000)	10	1000x PacBio Sequencing Program 2.0 (101-900-000) 10x Sequel II HiFi Read Cell (101-900-000) 10x Sequel II HiFi Read Cell (101-900-000) 10x Sequel II HiFi Read Cell (101-900-000)
De novo Assembly for Long-Read Phasing	Sequel II De Novo Phasing Long Reads-10 (101-900-000)	10	1000x PacBio Sequencing Program 2.0 (101-900-000) 10x Sequel II HiFi Read Cell (101-900-000) 10x Sequel II HiFi Read Cell (101-900-000) 10x Sequel II HiFi Read Cell (101-900-000)

[PacBio Application Consumable Bundle Purchasing Guide](#) (PG100-082620)

Purchasing Guide enables users to easily order required consumables needed to prepare a SMRTbell library to run a specific type of application on the Sequel II and IIe Systems.*



cDNA Synthesis & Multiplexed PCR Amplicon Generation
Amplify and asymmetrically barcode SARS-CoV-2 samples for multiplexing in a single library using barcoded M13 Primers

Low-Throughput and High-Throughput HiFiViral sample prep workflows available

Library Construction (SMRTbell Express TPK 2.0)

HiFi Sequencing (Sequel II and IIe Systems)

Recommend pooling up to 900 SARS-CoV-2 samples per SMRT Cell 8M

Data Analysis

Perform assembly or variant calling with HiFi data using Coronavirus (SARS-Cov-2) sequencing analysis (CoSA) tools on GitHub

* Note: For SARS-CoV-2 amplicon sequencing, users can choose the 'Amplicons for <3kb with Barcoded Primers' consumables bundle (101-901-000) and order other required kit products separately – refer to the [HiFiViral protocol](#) for specific details.

Detailed, End-to-End PacBio Protocol for SARS-CoV-2 Full-Viral Genome Sequencing



Highly accurate long reads offer critical advantages for surveillance



Customer Collaboration – PacBio HiFiViral for SARS-CoV-2 Workflow

High-Throughput Multiplexing 1.2 kb Amplicons for Full-Viral Genome Sequencing of SARS-CoV-2

This protocol is the result of a collaboration between LabCorp and Pacific Biosciences. We acknowledge the contributions of:

1. LabCorp
Brian Krueger, Michael Levandoski, Jonathan Williams, Brian Novvell, Scott Parker, Stanley Letovsky, Qian Zeng, Lax Iyer, Marcia Eisenberg
2. Pacific Biosciences
Primo Baybayan, George Yuan, Elizabeth Tseng, Jonas Kortach

This document describes a workflow for whole viral genome sequencing of SARS-CoV-2 samples on the Sequel[®] II and Sequel[®] Ite Systems using a targeted PCR approach. First-strand cDNA is synthesized using the SuperScript VIL0 cDNA Synthesis Kit (Thermo Fisher Scientific). The resulting cDNA is then PCR amplified, asymmetrically barcoded and subsequently pooled with other samples for SMRTbell[®] library construction and multiplex sequencing on the Sequel[®] II and Ite Systems.

If your throughput needs do not require multiplexing more than 48 samples, we recommend following the [Low-Throughput HiFiViral procedure](#) for full-viral genome sequencing of SARS-CoV-2. The [Low-Throughput workflow](#) is also recommended if you have little or no experience with high-throughput (e.g., >96-plex) sample preparation for multiplexed SMRT Sequencing. As you gain significant experience with PacBio's end-to-end workflow for multiplexed sample preparation, you may ramp up your throughput by using the High-Throughput procedure.

This procedure recommends amplification and sequencing of twenty-nine overlapping 1.2 kb amplicons that are tiled across the full 29.9 kb SARS-CoV-2 genome (Figure 1). The PCR primers were designed by Nikki E. Freed of the School of Natural and Computational Sciences, Massey University, Auckland, New Zealand.



Figure 1: 29 x 1.2 kb amplicons tiled across the SARS-CoV-2 genome.

This procedure requires two-rounds of amplification, first using M13-tailed target specific primers to tail the PCR products with a universal M13 sequence, followed by a second-round PCR using barcoded M13 primers. There are 32 Forward (F) and 32 Reverse (R) barcoded M13 primers available for barcoding and when used in different combinations allow multiplexing of up to 1024 samples.

[PacBio COVID-19 Sequencing Tools and Resources](#)



Easier Primer Balancing



Fewer Amplicon Dropouts



Catch All Variants



Flexible Batch Size



Cost Effective



PacBio HiFiViral for SARS-CoV-2 Workflow Overview

MULTIPLEXED SARS-CoV-2 AMPLICON LIBRARY PREPARATION PROCEDURE DESCRIPTION

- [Customer Collaboration – PacBio HiFiViral for SARS-CoV-2 Workflow](#) documents (PN [102-082-500](#) / [102-075-000](#)) describe a procedure for whole viral genome sequencing of asymmetrically barcoded SARS-CoV-2 samples on the Sequel II and IIe Systems using a 2-step targeted PCR amplicon approach.
- The procedure involves amplification and sequencing of twenty-nine overlapping 1.2 kb amplicons that are tiled across the full 29.9 kb SARS-CoV-2 genome
- This library preparation workflow is amenable to automation and was developed in collaboration with [Laboratory Corporation of America](#)
- Two versions of the HiFiViral protocol are available:
 1. **Low-throughput sample preparation workflow** (PN [102-082-500](#)) supports multiplexed library construction for up to 48 SARS-CoV-2 samples at a time
 2. **High-throughput sample preparation workflow** (PN [102-075-000](#)) supports multiplexed library construction for up to 900 SARS-CoV-2 samples at a time
- PacBio recommends pooling **up to 900 SARS-CoV-2 samples for sequencing on a single SMRT Cell 8M.**

The screenshot shows the title page of a document titled "Customer Collaboration – PacBio HiFiViral for SARS-CoV-2 Workflow". It includes the PacBio logo, a list of authors from LabCorp and Pacific Biosciences, and a detailed description of the workflow for whole viral genome sequencing. A diagram illustrates the 29.9 kb SARS-CoV-2 genome with 29 overlapping 1.2 kb amplicons. The document also mentions that PCR primers were designed by Nikki E. Freed of Massey University.

[PacBio COVID-19 Sequencing Tools and Resources](#)

RESEARCH FOCUS
MICROBIOLOGY AND INFECTIOUS DISEASE
 PacBio COVID-19 Sequencing Tools and Resources



PACBIO HIFIVIRAL FOR SARS-CoV-2 WORKFLOW OVERVIEW



cDNA Synthesis & Multiplexed PCR Amplicon Generation (~8 hours)

Amplify and asymmetrically barcode SARS-CoV-2 samples for multiplexing in a single library using PacBio-Barcoded M13 Primers
[Low-Throughput \($\leq 48\text{-plex}$\)](#) and [High-Throughput \(up to 900-plex\)](#) sample preparation workflows available



Library Construction (SMRTbell Express TPK 2.0) (~6 hours)




HiFi Sequencing (Sequel II and IIe Systems) (15-hour collection time)

Recommend pooling up to 900 SARS-CoV-2 samples per SMRT Cell 8M



Data Analysis

Perform assembly or variant calling with HiFi data using Coronavirus (SARS-Cov-2) sequencing analysis ([CoSA](#)) tools on GitHub




HIFI SEQUENCING HAS HIGHER SUCCESS IN DELIVERING COMPLETE VIRAL GENOMES

Easier Primer Balancing:

1.2 kb amplicons optimizes for even genome coverage by simplifying the PCR reaction while preserving performance with degraded samples



**Easier Primer
Balancing**




**Fewer Amplicon
Dropouts**

Fewer Amplicon Dropouts:

With fewer amplicons, primers anneal to less of the viral genome, lowering the chance that viral mutations will cause dropouts due to primer mismatches

Catch All Variants:



Longer, 99.9+% accurate HiFi reads detect both SNPs and indels up to 300 bp, allow phasing of variants





**Catch All
variants**

USING FEWER REQUIRED PRIMER SETS SIMPLIFIES PCR BALANCING

Unbalanced PCR reactions are a driver of uneven coverage, poor scale up, and missing regions



Balancing 96 or more primer sets




Balancing 29 primer sets for PacBio Sequencing

USING FEWER REQUIRED PRIMER SETS LEADS TO FEWER DROPOUTS


Less of the genome is covered by the primer sequences, making the assay more robust to coverage dropouts due to primer mismatches with new viral strains

29 amplicons




PacBio solution

98 amplicons




**Illumina / ONT ARTIC
solution**

237 amplicons




Ion Torrent solution

HIFI SARS-CoV-2 SEQUENCING ENABLES DETECTION OF LARGER VARIANTS



In collaboration with Labcorp

HIFI SARS-CoV-2 SEQUENCING ENABLES PHASING OF VARIANTS IN A SAMPLE CONTAINING MULTIPLE SUBTYPES



In collaboration with Labcorp

PACBIO'S HiFiViral SEQUENCING WORKFLOW PROVIDES THE RIGHT SCALE AT THE RIGHT COST PER SAMPLE

Flexible Batch Size:

The assay can be scaled to sequence from 48-900 samples per SMRT Cell on the Sequel IIe System

Cost Effective:

Per sample costs are well within average reimbursement rates for successfully recovered, complete sequences



**Flexible Batch
Size**



Cost Effective



Multiplexed SARS-CoV-2 Amplicon Library Preparation Using PacBio-Barcoded M13 Primers

PACBIO HiFiViral FOR SARS-CoV-2 WORKFLOW PROTOCOL DOCUMENT

- [Customer Collaboration – PacBio HiFiViral for SARS-CoV-2 Workflow](#) protocol describes a procedure for amplification and sequencing of twenty-nine overlapping 1.2 kb amplicons* that are tiled across the full 29.9 kb SARS-CoV-2 genome
 - **Low-Throughput procedure** (PN [102-082-500](#)) supports multiplexed library construction for sequencing up to 48 samples at a time
 - **High-Throughput procedure** (PN [102-075-000](#)) supports multiplexed library construction for sequencing up to 900 samples at a time
- Following first-strand cDNA synthesis, the resulting cDNA is then amplified using a 2-step PCR method, asymmetrically barcoded and subsequently pooled with other samples for SMRTbell library construction and multiplexed sequencing on the Sequel II or Ile System.

Protocol Document Contents

1. Link to [HiFiViral for SARS-CoV-2 Oligo Ordering Sheet](#) for ordering M13-tailed target-specific primers and barcoded M13 primers
2. Instructions for performing first-strand SARS-CoV-2 cDNA synthesis and amplification of cDNA products prior to SMRTbell library construction
3. Instructions for pooling amplified SARS-CoV-2 cDNA products and constructing SMRTbell libraries using SMRTbell Express Template Prep Kit 2.0
4. Sample setup guidance for preparing multiplexed SARS-CoV-2 amplicon SMRTbell libraries for sequencing on the Sequel II and Ile Systems

The thumbnail shows the cover page of the protocol document. It features the PacBio logo at the top left. The title is "Customer Collaboration – PacBio HiFiViral for SARS-CoV-2 Workflow High-Throughput Multiplexing 1.2 kb Amplicons for Full-Viral Genome Sequencing of SARS-CoV-2". Below the title, it lists the authors: LabCorp, Brian Khaviger, Michael Levandovski, Jonathan Williams, Brian Norvelt, Scott Pantzer, Sergey Letovsky, Qian Zeng, Lay Yee, Marcia Eisenberg, Pacific Biosciences, Pritam Bhatnagar, George Yuan, Elizabeth Tseng, Jonas Kortach. A short abstract follows, describing the workflow for whole viral genome sequencing of SARS-CoV-2 samples. A diagram shows a 29.9 kb genome with 29 overlapping 1.2 kb amplicons. The bottom of the page includes contact information for support and a footer with the document ID: 102-075-000 Version 03 (February 17, 2021).

[PacBio COVID-19 Sequencing Tools and Resources](#)

* 1.2 kb amplicon strategy for SARS-CoV-2 sequencing was adapted from Freed, N.E. et al. (2020) *Biology Methods & Protocols* 5(1):1 – 7

REQUIRED MATERIALS & EQUIPMENT

ITEM	WHERE USED	VENDOR	PART NUMBER
cDNA Synthesis			
SuperScript VILO cDNA Synthesis Kit	cDNA Preparation	Thermo Fisher Scientific	11754250
Oligo(dT) (100 μM)	cDNA Preparation	IDT	Custom*
PCR Reaction			
Target-Specific F/R Primers tailed with M13 Sequences (Customer-supplied)	PCR Amplification (1 st -Round)	Oligo Synthesis Company	N/A
F/R PacBio-Barcoded M13 Primers (Customer-supplied)	PCR Amplification (2 nd -Round)	Oligo Synthesis Company	N/A
Q5 Hot Start High-Fidelity DNA Polymerase	PCR Amplification (1 st -Round)	NEB	M0493L
dNTPS	PCR Amplification (1 st -Round)	NEB	N0447L
KAPA HiFi HotStart ReadyMix	PCR Amplification (2 nd -Round)	Roche	7958935001
Nuclease-Free Water	PCR Amplification	Any	Vendor-specific

* Oligo(dT) sequence is 5' TTT TTT TTT TTT GTC ATT CTC CTA AG 3' and HPLC purification is recommended

REQUIRED MATERIALS & EQUIPMENT (CONT.)

ITEM	WHERE USED	VENDOR	PART NUMBER
SMRTbell Library Construction			
SMRTbell Express Template Prep Kit 2.0	Library Preparation	PacBio	100-938-900
SMRTbell Enzyme Cleanup Kit	Library Preparation	PacBio	101-746-400
DynaMag-2 Magnet	Purification	Invitrogen	12321D
100% Ethanol, Molecular Biology Grade	Purification	Any	Vendor-specific
AMPure PB Beads	Purification	PacBio	100-265-900

PACBIO HiFiViral FOR SARS-CoV-2 WORKFLOW DETAILS




1. First-Strand cDNA Synthesis (~2 hrs)

- Perform 1st-strand cDNA synthesis using the SuperScript VILO cDNA Synthesis Kit.



2. cDNA Amplification and Pooling (~6 hrs)

- Prepare 1st and 2nd round PCR reactions to generate overlapping barcoded 1.2 kb amplicons tiled across the full 29.9 kb SARS-CoV-2 genome



- 1st-Round PCR

- Perform 1st-round PCR using SARS-CoV-2 specific primers tailed with M13 sequences
- Each sample requires that two multiplex PCR reactions using **Primer Pool 1 (15 primer pairs)** and **Primer Pool 2 (14 primer pairs)** be performed in parallel

- 2nd-Round PCR

- For the second-round PCR, PCR products from the 1st-round PCR reactions are re-amplified using barcoded M13 primers to generate **asymmetric barcoded amplicons**

- Pool Barcoded Amplicons

- Perform pooling of PacBio-barcoded amplicon products
- PacBio recommends pooling **up to 900 SARS-CoV-2 samples per SMRT Cell 8M**
- Purify single, pooled amplicon sample using AMPure PB beads

1

cDNA Preparation



1st- Round PCR with Target-Specific Primers Tailed with M13 Sequences

PCR Round 1



2nd-Round PCR with Asymmetric Barcoded M13 Primers (F/R)

PCR Round 2



Pool Barcoded Amplicons





3. SMRTbell Express TPK 2.0 Library Construction (~6 hrs)

- The amount of total pooled (barcoded) amplicon DNA required for SMRTbell library construction is 500 ng – 1000 ng.
- Typical library construction yield is $\geq 40\%$




4. HiFi Sequencing (15-hr collection time)

- Anneal sequencing primer (1 hour), bind polymerase (1 hour), and perform complex cleanup (0.5 hours)
- Use a 15-hour movie collection time for sequencing SARS-CoV-2 samples on the Sequel II or IIe System
- Generate HiFi reads with $>99.9\%$ accuracy



5. Data Analysis

- Perform assembly or variant calling with HiFi data using Coronavirus (SARS-Cov-2) sequencing analysis ([CoSA](#)) tools on GitHub



PREPARATION OF M13-TAILED SARS-CoV-2 PRIMERS FOR 1ST-ROUND PCR

- Forward and reverse SARS-CoV-2 PCR primers tailed with M13 sequences may be ordered from any oligo synthesis provider
- The [HiFiViral for SARS-CoV-2 Oligo Ordering Sheet](#) provides primer sequences (5' → 3') and recommended synthesis scale
 - Add a 5' blocker (e.g., 5AmMC6) to ensure that carryover amplicons from the 1st-Round PCR are not ligated to SMRTbell adapters during library construction. [Desalted primers are sufficient for 1st-Round PCR amplification.]
- Before use, dilute oligos with nuclease-free buffer (10mM Tris-HCl pH 7.5) to 100 μM stock concentration.
- Prepare **100 μM** stocks of **Primer Pool 1 (odd primer pairs)** and **Primer Pool 2 (even primer pairs)** by adding the appropriate volumes of F and R primers as shown in the tables below

Primer Pool 1 (250 μL): Odd Primer Pairs (15 Oligo Pairs (F/R); 30 Oligos Total)

Primer Pair	Forward Primer to Add (μL)	Reverse Primer to Add (μL)
SARSCoV_1200_1	10.0	10.0
SARSCoV_1200_3	10.0	10.0
SARSCoV_1200_5	5.0	5.0
SARSCoV_1200_7	10.0	10.0
SARSCoV_1200_9	10.0	10.0
SARSCoV_1200_11	5.0	5.0
SARSCoV_1200_13	10.0	10.0
SARSCoV_1200_15	10.0	10.0
SARSCoV_1200_17	5.0	5.0
SARSCoV_1200_19	10.0	10.0
SARSCoV_1200_21	5.0	5.0
SARSCoV_1200_23	5.0	5.0
SARSCoV_1200_25	10.0	10.0
SARSCoV_1200_27	10.0	10.0
SARSCoV_1200_29	10.0	10.0

Primer Pool 2 (200 μL): Even Primer Pairs (14 Oligo Pairs (F/R); 28 Oligos Total)

Primer Pair	Forward Primer to Add (μL)	Reverse Primer to Add (μL)
SARSCoV_1200_2	5.0	5.0
SARSCoV_1200_4	5.0	5.0
SARSCoV_1200_6	5.0	5.0
SARSCoV_1200_8	10.0	10.0
SARSCoV_1200_10	10.0	10.0
SARSCoV_1200_12	5.0	5.0
SARSCoV_1200_14	10.0	10.0
SARSCoV_1200_16	5.0	5.0
SARSCoV_1200_18	5.0	5.0
SARSCoV_1200_20	5.0	5.0
SARSCoV_1200_22	10.0	10.0
SARSCoV_1200_24	10.0	10.0
SARSCoV_1200_26	5.0	5.0
SARSCoV_1200_28	10.0	10.0

PREPARATION OF PACBIO-BARCODED M13 PRIMERS FOR 2ND-ROUND PCR

- 32 forward and 32 reverse PacBio-barcoded M13 primer sequences are available for the 2nd-round PCR step to create asymmetrically barcoded SARS-CoV-2 amplicons
 - PacBio recommends multiplexing **up to 900 SARS-CoV-2 samples per SMRT Cell 8M** for HiFi sequencing
- The [HiFiViral for SARS-CoV-2 Oligo Ordering Sheet](#) provides primer sequences (5' → 3') and recommended synthesis scale
 - HPLC-purified primers are recommended for 2nd-Round PCR amplification
- Resuspend PacBio-barcoded M13 primers with nuclease-free buffer (10mM Tris-HCl pH 7.5) to a concentration of 3.0 μM and aliquot into a 96-well plate as shown in the example layout below

EXAMPLE PLATE LAYOUT FOR RESUSPENDING 32 FORWARD AND 32 REVERSE PACBIO-BARCODED M13 PRIMER OLIGOS.

	1	2	3	4	5	6	7	8	9	10	11	12
A	FWD_1001	FWD_1009	FWD_1017	FWD_1025	x	x	REV_1049	REV_1057	REV_1065	REV_1073	x	x
B	FWD_1002	FWD_1010	FWD_1018	FWD_1026	x	x	REV_1050	REV_1058	REV_1066	REV_1074	x	x
C	FWD_1003	FWD_1011	FWD_1019	FWD_1027	x	x	REV_1051	REV_1059	REV_1067	REV_1075	x	x
D	FWD_1004	FWD_1012	FWD_1020	FWD_1028	x	x	REV_1052	REV_1060	REV_1068	REV_1076	x	x
E	FWD_1005	FWD_1013	FWD_1021	FWD_1029	x	x	REV_1053	REV_1061	REV_1069	REV_1077	x	x
F	FWD_1006	FWD_1014	FWD_1022	FWD_1030	x	x	REV_1054	REV_1062	REV_1070	REV_1078	x	x
G	FWD_1007	FWD_1015	FWD_1023	FWD_1031	x	x	REV_1055	REV_1063	REV_1071	REV_1079	x	x
H	FWD_1008	FWD_1016	FWD_1024	FWD_1032	x	x	REV_1056	REV_1064	REV_1072	REV_1082	x	x

Columns 1-4 are M13 forward primers tailed with PacBio barcode 1001 to barcode 1032.

Columns 7-10 are M13 reverse primers tailed with PacBio barcode 1049 to barcode 1079 and barcode 1082.

EXAMPLE ASYMMETRIC BARCODE PLATE MAPS FOR HIGH-THROUGHPUT MULTIPLEXING OF SARS-CoV-2 SAMPLES

- See **Appendix 1** in the high-throughput procedure (PN [102-075-000](#)) for example asymmetric barcode plate maps for processing up to 10 x 96-well plates containing SARS-Cov-2 samples at a time

Barcode Plate 1

	1	2	3	4	5	6	7	8	9	10	11	12
A	1001	1009	1017	1025	1001	1009	1017	1025	1001	1009	1017	1025
B	1008	1049	1041	1068	1052	1050	1058	1052	1050	1058	1052	1050
C	1003	1011	1019	1027	1003	1011	1019	1027	1003	1011	1019	1027
D	1004	1042	1030	1028	1004	1042	1030	1028	1004	1042	1030	1028
E	1006	1044	1040	1048	1006	1044	1040	1048	1006	1044	1040	1048
F	1005	1013	1021	1029	1005	1013	1021	1029	1005	1013	1021	1029
G	1002	1040	1042	1049	1002	1040	1042	1049	1002	1040	1042	1049
H	1007	1015	1023	1031	1007	1015	1023	1031	1007	1015	1023	1031
I	1009	1048	1041	1049	1009	1048	1041	1049	1009	1048	1041	1049
J	1008	1049	1049	1049	1008	1049	1049	1049	1008	1049	1049	1049

Barcode Plate 2

	1	2	3	4	5	6	7	8	9	10	11	12
A	1001	1009	1017	1025	1001	1009	1017	1025	1001	1009	1017	1025
B	1002	1010	1018	1026	1002	1010	1018	1026	1002	1010	1018	1026
C	1003	1011	1019	1027	1003	1011	1019	1027	1003	1011	1019	1027
D	1004	1012	1020	1028	1004	1012	1020	1028	1004	1012	1020	1028
E	1005	1013	1021	1029	1005	1013	1021	1029	1005	1013	1021	1029
F	1006	1014	1022	1030	1006	1014	1022	1030	1006	1014	1022	1030
G	1007	1015	1023	1031	1007	1015	1023	1031	1007	1015	1023	1031
H	1008	1016	1024	1032	1008	1016	1024	1032	1008	1016	1024	1032
I	1009	1017	1025	1033	1009	1017	1025	1033	1009	1017	1025	1033
J	1008	1016	1024	1032	1008	1016	1024	1032	1008	1016	1024	1032

Barcode Plate 3

	1	2	3	4	5	6	7	8	9	10	11	12
A	1001	1009	1017	1025	1001	1009	1017	1025	1001	1009	1017	1025
B	1002	1010	1018	1026	1002	1010	1018	1026	1002	1010	1018	1026
C	1003	1011	1019	1027	1003	1011	1019	1027	1003	1011	1019	1027
D	1004	1012	1020	1028	1004	1012	1020	1028	1004	1012	1020	1028
E	1005	1013	1021	1029	1005	1013	1021	1029	1005	1013	1021	1029
F	1006	1014	1022	1030	1006	1014	1022	1030	1006	1014	1022	1030
G	1007	1015	1023	1031	1007	1015	1023	1031	1007	1015	1023	1031
H	1008	1016	1024	1032	1008	1016	1024	1032	1008	1016	1024	1032
I	1009	1017	1025	1033	1009	1017	1025	1033	1009	1017	1025	1033
J	1008	1016	1024	1032	1008	1016	1024	1032	1008	1016	1024	1032

Barcode Plate 4

	1	2	3	4	5	6	7	8	9	10	11	12
A	1001	1009	1017	1025	1001	1009	1017	1025	1001	1009	1017	1025
B	1002	1010	1018	1026	1002	1010	1018	1026	1002	1010	1018	1026
C	1003	1011	1019	1027	1003	1011	1019	1027	1003	1011	1019	1027
D	1004	1012	1020	1028	1004	1012	1020	1028	1004	1012	1020	1028
E	1005	1013	1021	1029	1005	1013	1021	1029	1005	1013	1021	1029
F	1006	1014	1022	1030	1006	1014	1022	1030	1006	1014	1022	1030
G	1007	1015	1023	1031	1007	1015	1023	1031	1007	1015	1023	1031
H	1008	1016	1024	1032	1008	1016	1024	1032	1008	1016	1024	1032
I	1009	1017	1025	1033	1009	1017	1025	1033	1009	1017	1025	1033
J	1008	1016	1024	1032	1008	1016	1024	1032	1008	1016	1024	1032

Barcode Plate 5

	1	2	3	4	5	6	7	8	9	10	11	12
A	1001	1009	1017	1025	1001	1009	1017	1025	1001	1009	1017	1025
B	1002	1010	1018	1026	1002	1010	1018	1026	1002	1010	1018	1026
C	1003	1011	1019	1027	1003	1011	1019	1027	1003	1011	1019	1027
D	1004	1012	1020	1028	1004	1012	1020	1028	1004	1012	1020	1028
E	1005	1013	1021	1029	1005	1013	1021	1029	1005	1013	1021	1029
F	1006	1014	1022	1030	1006	1014	1022	1030	1006	1014	1022	1030
G	1007	1015	1023	1031	1007	1015	1023	1031	1007	1015	1023	1031
H	1008	1016	1024	1032	1008	1016	1024	1032	1008	1016	1024	1032
I	1009	1017	1025	1033	1009	1017	1025	1033	1009	1017	1025	1033
J	1008	1016	1024	1032	1008	1016	1024	1032	1008	1016	1024	1032

Barcode Plate 6


	1	2	3	4	5	6	7	8	9	10	11	12
A	1001	1009	1017	1025	1001	1009	1017	1025	1001	1009	1017	1025
B	1002	1010	1018	1026	1002	1010	1018	1026	1002	1010	1018	1026
C	1003	1011	1019	1027	1003	1011	1019	1027	1003	1011	1019	1027
D	1004	1012	1020	1028	1004	1012	1020	1028	1004	1012	1020	1028
E	1005	1013	1021	1029	1005	1013	1021	1029	1005	1013	1021	1029
F	1006	1014	1022	1030	1006	1014	1022	1030	1006	1014	1022	1030
G	1007	1015	1023	1031	1007	1015	1023	1031	1007	1015	1023	1031
H	1008	1016	1024	1032	1008	1016	1024	1032	1008	1016	1024	1032
I	1009	1017	1025	1033	1009	1017	1025	1033	1009	1017	1025	1033
J	1008	1016	1024	1032	1008	1016	1024	1032	1008	1016	1024	1032

Barcode Plate 10


	1	2	3	4	5	6	7	8	9	10	11	12
A	1001-1076	1009-1076	1017-1076	1025-1076	1001-1077	1009-1077	1017-1077	1025-1077	1001-1078	1009-1078	1017-1078	1025-1078
B	1002-1076	1010-1076	1018-1076	1026-1076	1002-1077	1010-1077	1018-1077	1026-1077	1002-1078	1010-1078	1018-1078	1026-1078
C	1003-1076	1011-1076	1019-1076	1027-1076	1003-1077	1011-1077	1019-1077	1027-1077	1003-1078	1011-1078	1019-1078	1027-1078
D	1004-1076	1012-1076	1020-1076	1028-1076	1004-1077	1012-1077	1020-1077	1028-1077	1004-1078	1012-1078	1020-1078	1028-1078
E	1005-1076	1013-1076	1021-1076	1029-1076	1005-1077	1013-1077	1021-1077	1029-1077	1005-1078	1013-1078	1021-1078	1029-1078
F	1006-1076	1014-1076	1022-1076	1030-1076	1006-1077	1014-1077	1022-1077	1030-1077	1006-1078	1014-1078	1022-1078	1030-1078
G	1007-1076	1015-1076	1023-1076	1031-1076	1007-1077	1015-1077	1023-1077	1031-1077	1007-1078	1015-1078	1023-1078	1031-1078
H	1008-1076	1016-1076	1024-1076	1032-1076	1008-1077	1016-1077	1024-1077	1032-1077	1008-1078	1016-1078	1024-1078	1032-1078

EXAMPLE **HIGH-THROUGHPUT*** SAMPLE PREPARATION WORKFLOW FOR FULL-VIRAL GENOME SEQUENCING OF SARS-CoV-2 SAMPLES

1. Perform one 1st-strand cDNA synthesis reaction per sample
2. Perform 1st-round PCR using two separate, parallel PCR reactions per sample with SARS-CoV-2 PCR primers tailed with M13 sequences
 - **Primer Pool 1** contains 15 **odd** primer pairs
 - **Primer Pool 2** contains 14 **even** primer pairs




3. Perform 2nd-round PCR using two separate, parallel PCR reactions per sample with PacBio-Barcoded M13 Primers
 - See **Appendix 1** in the high-throughput procedure for example asymmetric barcode plate maps for processing up to 10 x 96-well plates containing SARS-CoV-2 samples at a time
4. Combine 2nd-round PCR products (1- μ L aliquots) from each well of R2 Pool 1 and R2 Pool 2 into a single 2 mL LoBind tube and perform AMPure PB bead purification prior to proceeding with SMRTbell library construction




* See *HiFiViral Workflow for High-Throughput Multiplexing of 1.2 kb Amplicons for Full-Viral Genome Sequencing of SARS-CoV-2* (PN [102-075-000](#))

EXAMPLE **HIGH-THROUGHPUT** SAMPLE PREPARATION WORKFLOW FOR FULL-VIRAL GENOME SEQUENCING OF SARS-CoV-2 SAMPLES (CONT.)




EXAMPLE LOW-THROUGHPUT* SAMPLE PREPARATION WORKFLOW FOR FULL-VIRAL GENOME SEQUENCING OF SARS-CoV-2 SAMPLES

1. Perform one 1st-strand cDNA synthesis reaction per sample
2. Perform 1st-round PCR using two separate, parallel PCR reactions per sample with SARS-CoV-2 PCR primers tailed with M13 sequences
 - **Primer Pool 1** contains 15 **odd** primer pairs
 - **Primer Pool 2** contains 14 **even** primer pairs



3. Perform 2nd-round PCR using two separate, parallel PCR reactions per sample with PacBio-Barcoded M13 Primers
 - See **Appendix 1** in the low-throughput procedure for an example asymmetric barcode plate map for processing up to 48 SARS-Cov-2 samples in a single 96-well plate
4. Normalize R2 Pool 1 and R2 Pool 2 PCR products to 10 nM to obtain more balanced coverage.
5. Combine normalized 2nd-round PCR products (2- μ L aliquots from each well) into a single 2 mL LoBind tube and perform AMPure PB bead purification prior to proceeding with SMRTbell library construction



* See *HiFiViral Workflow for Low-Throughput Multiplexing of 1.2 kb Amplicons for Full-Viral Genome Sequencing of SARS-CoV-2* (PN [102-082-500](#))



Multiplexed SARS-CoV-2 Amplicon Library Sequencing Workflow Recommendations

SAMPLE SETUP RECOMMENDATIONS FOR SARS-CoV-2 AMPLICON LIBRARIES – SEQUEL II AND IIe SYSTEMS

- Follow **SMRT Link Sample Setup** instructions using the recommendations provided in the tables below for sequencing SARS-CoV-2 amplicon samples.
 - For **SMRT Link v10.0** (or higher): Select '**Amplicons <3 kb**' from the **Application** field drop-down menu in the SMRT Link Sample Setup and SMRT Link Run Design user interface and enter in the values shown in the tables below

SAMPLE SETUP CONDITIONS	SEQUEL II AND IIe SYSTEMS
Sequencing Primer	Sequencing Primer v4
Primer to Template Ratio	20:1
Polymerase to Template Ratio	10:1
Binding Kit	Sequel II Binding Kit 2.1
Binding Time	1 hour
Complex Cleanup Method	AMPure PB Beads
AMPure PB Bead Cleanup Anticipated Yield	35%


RUN DESIGN CONDITIONS	SEQUEL II AND IIe SYSTEMS
Sequencing Kit	Sequel II Sequencing Plate 2.0
Recommended On-Plate Loading Concentration	100 pM – 160 pM
Movie Collection Time	15 hours
Pre-extension Time	1 hour



Multiplexed SARS-CoV-2 Amplicon Data Analysis Recommendations

SARS-CoV-2 SEQUENCING DATA ANALYSIS TOOLS

Coronavirus (SARS-CoV-2) Sequencing Analysis (CoSA) is a set of Python and R scripts available on GitHub for analyzing SARS-CoV-2 sequences from PacBio HiFi / CCS data



PacificBiosciences / CoSA

Code Issues Pull requests Actions Projects Wiki Security Insights

master 2 branches 4 tags

Go to file Code

Magdoll v8.3.0 f6f3fb2 3 hours ago 93 commits

R	v2.0.0 submit	10 months ago
beta	v6.1.0	5 days ago
build	v8.2.0	15 hours ago
cosa	v6.1.0	5 days ago
data	adding freed protocol primer	5 days ago
examples	v6.1.0	5 days ago
vcf	v8.3.0	3 hours ago
LICENSE.md	commit cleanup for v4.0.0	11 days ago
README.md	v8.3.0	3 hours ago
setup.py	v8.3.0	3 hours ago

README.md

CoSA

Coronavirus (SARS-Cov-2) sequencing analysis

Last Updated: 02.22.2020 (v8.3.0)

Stable version of CoSA under PacBio's GitHub

About

SARS-CoV-2 analysis using PacBio long reads

sequencing pacbio sars-cov-2

Readme

View license

Releases

4 tags

Packages

No packages published

Contributors 3

- Magdoll Elizabeth Tseng
- jrharting John Harting
- natechols Nat Echols

Languages

Python 90.7% R 5.5%

Join the **COVID-19 Google Group** for CoSA updates

PREREQUISITES FOR SARS-CoV-2 VARIANT CALLING ANALYSIS USING PACBIO HIFI DATA

Python Requirements

- Python 3.7+
- [BioPython](#)
- [PyVCF](#)
- [mappy](#) and [panda](#) (if using [pbaa](#))

We recommend looking into [Anaconda](#) for easily managing Python packages

PacBio Tool Prerequisites

- [ccs](#)
- [lima](#)
- optional - [pbaa](#) if using *pbaa* for variant calling.


These PacBio tools can be obtained either through [SMRT Link](#) or [pbbioconda](#)

Other Prerequisites

- [CoSA](#)
- [samtools](#)
- [bamtools](#)
- optional - [bcftools](#), if using *bcftools* for variant calling.

GENERAL WORKFLOW FOR SARS-CoV-2 VARIANT CALLING ANALYSIS

1. Generate CCS data*
2. De-multiplex different barcoded samples (If dataset is multiplexed)
3. Trim PCR amplicon primers
4. Perform variant calling using:
 - [bcftools](#); or
 - [DeepVariant](#); or
 - [pbaa](#)
5. Generate consensus sequences
 - Use [VCFCons.py](#)
6. Assign lineages
 - Use [Pangolin](#) or [Nextclade](#)



Tool Choices	Output Format
--------------	---------------

<i>lima</i>	BAM
-------------	-----

<i>bcftools</i> <i>DeepVariant</i> <i>pbaa</i>	VCF
--	-----

<i>VCFCons</i> <i>(part of CoSA)</i>	FASTA
---	-------

* CCS analysis can be run using command line tools or by using the "Circular Consensus Sequencing (CCS)" application in [SMRT Link](#). For Sequel IIe Systems, CCS analysis can also be performed on-instrument.

EXAMPLE SARS-CoV-2 VARIANT CALLING ANALYSIS OUTPUT

```

CLC
NC_045512v2 8317 . A AC 200.00 . . GT:CLCAD2:DP 0/1:165,28:193
NC_045512v2 9867 . T C 200.00 . . GT:CLCAD2:DP 1/1:0,152:153
NC_045512v2 10450 . C T 200.00 . . GT:CLCAD2:DP 1/1:0,423:434
NC_045512v2 11287 . GTCTGGTTTT G 200.00 . . GT:CLCAD2:DP 1/1:0,424:435

bcftools
NC_045512v2 9867 . T C 228 . DP=154;VDB=0;SGB=-0.693147;RPB=1;MQB=1;MQSB=1;BQB=1;MQ0F=0;AC=2;AN=2;DP4=0,1,56,94;MQ=60 GT:PL 1/1:255,255,0
NC_045512v2 10450 . C T 228 . DP=99;VDB=0;SGB=-0.693147;RPB=0.999369;MQB=1;BQB=0.805608;MQ0F=0;AC=2;AN=2;DP4=4,0,95,0;MQ=60 GT:PL 1/1:255,108,0
NC_045512v2 11287 . GTCTGGTTTT G 228 . INDEL;IDV=94;IMF=0.949495;DP=99;VDB=0;SGB=-0.693147;MQ0F=0;AC=2;AN=2;DP4=4,0,95,0;MQ=60 GT:PL 1/1:255,156,0


DeepVariant
NC_045512v2 8317 . A AC 0 RefCall . GT:GQ:DP:AD:VAF:PL 0/0:47:200:168,30:0.15:0,46,64
NC_045512v2 9867 . T C 83.4 PASS . GT:GQ:DP:AD:VAF:PL 1/1:61:153:1,152:0.993464:83,61,0
NC_045512v2 10450 . C T 50.2 PASS . GT:GQ:DP:AD:VAF:PL 1/1:33:435:8,423:0.972414:50,33,0
NC_045512v2 11287 . GTCTGGTTTT G 41.2 PASS . GT:GQ:DP:AD:VAF:PL 1/1:27:431:7,424:0.983759:41,27,0

pbaa
NC_045512v2 9867 . T C . PASS NS=1;AF=1 GT:DP:AQ:AD:VAF:TG:HP:DV:CH 1:154:63.0495:153:1:NC_045512_9351_10400:0:-1:-1
NC_045512v2 10450 . C T . PASS NS=1;AF=1 GT:DP:AQ:AD:VAF:TG:HP:DV:CH 1:435:62.9768:424:1:NC_045512_10401_11420:0:-1:-1
NC_045512v2 11287 . GTCTGGTTTT G . PASS NS=1;AF=1 GT:DP:AQ:AD:VAF:TG:HP:DV:CH 1:435:62.9768:424:1:NC_045512_10401_11420:0:-1:-1
  
```

↓
VCFCons.py
 ↓

<p>Filtered variant calls</p> <pre> NC_045512v2 9867 . T C . PASS NC_045512v2 10450 . C T . PASS NC_045512v2 11287 . GTCTGGTTTT G . PASS </pre>	<p>Filtered consensus sequence</p> <pre> NC_045512v2T...C...GTCTGGTTTTAAGCT...TTATGG... SampleC...T...GAAGCT.....NNNNNN... </pre>
--	--

EXAMPLE SARS-CoV-2 VARIANT CALLING ANALYSIS VISUALIZATION



- Example visualization of SARS-Cov-2 variant calling analysis results using Integrative Genomics Viewer ([IGV](#))



Technical Documentation & Applications Support Resources

TECHNICAL RESOURCES FOR SARS-CoV-2 LIBRARY PREPARATION, SEQUENCING & DATA ANALYSIS

Visit PacBio's [COVID-19 Sequencing Tools and Resources Website](#) for HiFiViral for SARS-CoV-2 Workflow Updates and Other Resources

Sample Preparation Literature

- Customer Collaboration – PacBio HiFiViral for SARS-CoV-2 Workflow Protocol References
 - [Low-Throughput Procedure for Multiplexing 1.2 kb Amplicons for Full-Viral Genome Sequencing of SARS-CoV-2](#) (PN 102-082-500)
 - [High-Throughput Procedure for Multiplexing 1.2 kb Amplicons for Full-Viral Genome Sequencing of SARS-CoV-2](#) (PN 102-075-000)
- [Quick Reference Card – Loading and Pre-extension Recommendations for the Sequel II/IIe Systems](#) PN 101-769-100)
- [Overview – Sequel Systems Application Options and Sequencing Recommendations](#) (PN 101-851-300)
- [Application Consumable Bundles Purchasing Guide](#) (PN PG100-051320)
- [PacBio HiFiViral Workflow Overview: Multiplexed Amplicon Library Preparation for Full-Viral Genome Sequencing of SARS-CoV-2](#) (PN 102-084-800)

TECHNICAL RESOURCES FOR SARS-CoV-2 LIBRARY PREPARATION, SEQUENCING & DATA ANALYSIS (CONT.)

Data Analysis Resources

- Coronavirus (SARS-Cov-2) Sequencing Analysis ([CoSA](#)) GitHub Page
 - CoSA is a set of Python and R scripts for analyzing SARS-CoV-2 sequences from PacBio HiFi / CCS data.
- Join the [COVID19 Google Group](#) to stay up-to-date on bioinformatics analysis recommendations and [CoSA](#) changes

Videos & Webinars

- ASM 2020 Presentation (2020): Geographic and Temporal Mapping of the SARS-CoV-2 Pandemic in the United States [[Webinar Recording](#)]
- PacBio LabRoots Webinar (2020): Opportunities for using PacBio Long-read Sequencing for COVID-19 Research [[Webinar Recording](#)]
- PacBio Webinar (2020): Understanding SARS-CoV-2 and Host Immune Response to COVID-19 with PacBio Sequencing [[Webinar Recording](#)]



Appendix: RNA Isolation Kit Options for Full-Viral Genome Sequencing of SARS-CoV-2

RNA SAMPLE EXTRACTION KIT OPTIONS FOR FULL-VIRAL GENOME SEQUENCING OF SARS-CoV-2

Note: The products below have **not** been tested or validated by PacBio but are listed here as examples of third-party kits used by other PacBio customers for isolating SARS-CoV-2 RNA samples for multiplexed SMRTbell amplicon library preparation

VENDOR	RNA ISOLATION KIT PRODUCT	AUTOMATION PLATFORM
Thermo Fisher Scientific	MagMAX Viral and Pathogen Nucleic Acid Isolation Kit (Link)	KingFisher Flex System
Roche Molecular Systems	MagNA Pure 96 DNA and Viral NA Small Volume Kit (Link)	Roche MagNA Pure-96 (MP6)



PACBIO®

www.pacb.com

For Research Use Only. Not for use in diagnostic procedures. © Copyright 2021 by Pacific Biosciences of California, Inc. All rights reserved. Pacific Biosciences, the Pacific Biosciences logo, PacBio, SMRT, SMRTbell, Iso-Seq, and Sequel are trademarks of Pacific Biosciences. Pacific Biosciences does not sell a kit for carrying out the overall No-Amp Targeted Sequencing method. Use of these No-Amp methods may require rights to third-party owned intellectual property. FEMTO Pulse and Fragment Analyzer are trademarks of Agilent Technologies Inc.

All other trademarks are the sole property of their respective owners.