

Amplifying genomic DNA for SMRTbell library preparation and HiFi sequencing


Procedure & checklist

Overview

This procedure describes the workflow for constructing SMRTbell[®] libraries with PCR amplification from as little as 1 ng of starting genomic DNA. Sequencing libraries prepared with this procedure are compatible with PacBio[®] Revio[®] and Vega[™] systems.

Protocol overview		
Application	Genome sequencing	This protocol is intended for samples where there is insufficient DNA for PCR-free library prep, where difficult-to-remove sequencing inhibitors may be present, or where there is extensive DNA damage present
Genomic DNA Input	1–50 ng per sample	Recommended starting double-stranded DNA input amount for protocol
Target DNA fragment size	7–11 kb	Average fragment size of sheared DNA to balance read length with amplification efficiency
PCR polymerase	KOD Xtreme Hot Start DNA polymerase	Minimizes GC coverage bias with efficient amplification
Amplification adapter	Twist Universal Adapter System	Common adapter for NGS that enables unique dual-indexed (UDI) barcoding prior to SMRTbell library prep
SMRTbell library prep	SMRTbell prep kit 3.0	Supports the ligation of amplification adapters and the SMRTbell library preparation of amplified DNA for up to 24 samples
	Revio (SPRQ[™] chemistry)	Revio (non-SPRQ chemistry)/Vega
Required mass of amplified DNA into library prep per SMRT [®] Cell	150 ng	600 ng
	Time	
Protocol time from gDNA shearing to SMRTbell library	≤7 hr	
Protocol time for ABC workflow (polymerase binding)	1 hr	

Workflow overview



Required materials and equipment

DNA quantification	
Qubit fluorometer	ThermoFisher Scientific Q33238
Qubit 1X dsDNA HS assay kit	ThermoFisher Scientific Q33230
DNA shearing	
Megaruptor 3	Diagenode B06010003
Megaruptor 3 shearing kit	Diagenode E07010003
<i>or</i>	
g-TUBE	Covaris 520079 or 520104
PCR amplification and SMRTbell® library preparation	
KOD Xtreme™ Hot Start DNA polymerase (200U) (outside of Japan) KOD FX (200U) (Japan)	Millipore Sigma 71975-3* TOYOBO KFX-101
Twist Universal Adapter System	Twist 101307-101311
SMRTbell® prep kit 3.0, includes: <ul style="list-style-type: none"> • SMRTbell® prep kit 3.0 • SMRTbell® cleanup beads • Low TE buffer 	PacBio® 102-182-700
AMPure® PB bead size selection kit	PacBio® 102-182-500
Revio® SPRQ™ polymerase kit <i>or</i>	PacBio® 103-496-900
Vega™ polymerase kit <i>or</i>	PacBio® 103-426-500
Revio® polymerase kit (non-SPRQ™)**	PacBio® 102-739-100
*Manufactured by Toyobo, Japan and distributed by Merck KGaA globally.	
**Procedure for Revio polymerase kit (non-SPRQ) can be found in SMRT® Link Sample Setup	
Other common lab supplies	
100% ethanol, molecular biology	Any major lab supplier (MLS)
0.2 mL 8-tube strips	USA Scientific TempAssure 1402-4708
Magnetic separation rack compatible with 0.2 mL 8-tube strips	V&P Scientific VP 772F4-1
Thermocycler compatible with 100 µL and 8-tube strips	Any MLS
Microcentrifuge	Any MLS

Before you begin

The following are best practices for preparing SMRTbell libraries using PCR for HiFi sequencing on Revio or Vega systems. Please read carefully prior to beginning the procedure.

Application

This procedure is intended for, but not limited to, the following sample types:

- Those that do not have sufficient DNA mass for PCR-free sequencing (e.g., a single, small insect)
- Samples contaminated with sequencing inhibitors that are challenging to remove (e.g., snails or marine organisms)
- Preserved specimens collected in the field (e.g., ethanol preserved tissue)
- Chromosomal conformation capture (3C) samples (e.g., CiFi application)
- DNA from FFPE tissue (please note that DNA quantity and quality is highly variable from FFPE)

The mean library size is expected to be less than 10 kb due to the limitations of PCR. Please note that methylation signals will not be retained with this application.

Genomic DNA (gDNA) QC and input amount recommendations

Sequencing yield expectations should be set in accordance with the genomic DNA quality. Samples with higher amounts of degraded or fragmented DNA will have lower average insert sizes due to preferential amplification of the shorter fragments. To maximize HiFi base yield and average read length, it is recommended to start with HMW DNA. If longer HiFi read lengths are required for more degraded samples, pre- or post-library size selection methods may be used.

gDNA quality QC

The Agilent Femto Pulse system is highly recommended for accurately measuring the distribution gDNA mass. Please see the PacBio [Technical note](#) for more details.

Recommended guidelines for evaluating gDNA quality for this protocol:

- Use the Femto Pulse gDNA 165 kb analysis kit (Agilent FP-1002-0275)
- Dilute aliquot of the sample to 250 pg/ μ L

gDNA input mass

The supported genomic DNA input mass for this protocol ranges from 1 to 50 ng. The number of PCR cycles required will depend on the input mass.

Size selection of SMRTbell library

This protocol uses the diluted 3.1X (35% v/v) AMPure PB size selection and cleanup method to progressively deplete DNA below 5 kb. For libraries with an average size less than 5 kb, a standard 1X SMRTbell bead cleanup is recommended (e.g., FFPE samples)

Gel size selection approaches can be used to enrich for larger insert sizes when there is sufficient library available. Please see this [Technical Note](#) on gel cassette size selection.

Multiplexing

Please visit our [multiplexing webpage](#) for the demultiplexing fasta file to import into SMRT® Link versions 25.1 and above. We recommend using this file for demultiplexing and adapter trimming.

Reagent and sample handling

Room temperature is defined as any temperature in the range of **18–25°C** for this protocol.

SMRTbell cleanup beads and AMPure PB beads

Bring SMRTbell cleanup beads and AMPure PB beads to room temperature prior to use.

Vortex *immediately before* any addition to sample. Failure to do this will result in low recovery.

Pipette-mix all bead binding and elution steps until beads are distributed evenly in solution.

SMRTbell prep kit 3.0

Thaw the Repair buffer, Nuclease buffer, and SMRTbell adapter at room temperature. Mix reagent buffers with a brief vortex prior to use.

Thaw Elution buffer at room temperature. Elution buffer can be left at room temperature long term after thawing.

Quick-spin all reagents in microcentrifuge to collect liquid at bottom prior to use.

Keep all temperature-sensitive reagents on ice.

Temperature-sensitive reagents	
Reagent	Tube color
End repair mix	Blue
DNA repair mix	Green
SMRTbell adapter	Orange
Ligation mix	Yellow
Ligation enhancer	Red
Nuclease mix	Light green

Samples can be stored at 4°C at all safe stopping points listed in the protocol.

Anneal, bind, and cleanup using the Revio or Vega Polymerase Kit

Thaw the following reagents at room temperature:

Component	Tube color
Annealing buffer	Light blue
Standard sequencing primer	Light green
Polymerase buffer	Yellow
Loading buffer	Green
Dilution buffer	Blue

Once thawed, reaction buffers and sequencing primer may be stored on a cold block, at 4°C, or on-ice. The Loading buffer should be left at room-temperature.

Note: The Loading buffer is light sensitive and should be protected from light when not in use.

Keep the following reagents on a cold block or ice:

- Sequencing polymerase
- Sequencing control

Bring the following reagents to room temperature 30 minutes prior to use:

- Loading buffer

Third party reagents

Thaw 2x Xtreme buffer, dNTP, Twist universal adapter and UDI primers at room temperature.

Safety precautions

Refer to the Safety Data Sheet (SDS) for information on reagent hazards and protocols for safe handling, use, storage, and disposal.

Workflow steps

1. DNA shearing

We recommend shearing to a size range of 7–11 kb using either Diagenode Megaruptor 3 system (option 1) or Covaris g-TUBE (option 2). An example of the recommended shearing size distribution is shown in Figure 1.

Option 1: DNA shearing with the Megaruptor 3 system

✓	Step	Instructions				
	1.1	Bring 1–50 ng DNA up to a final volume of 65 µL with low TE buffer .				
		Shear DNA on the Megaruptor 3 system with the following parameters:				
		<table border="1"> <thead> <tr> <th>Shear speed</th> <th>Target insert length</th> </tr> </thead> <tbody> <tr> <td>59</td> <td>10 kb</td> </tr> </tbody> </table>	Shear speed	Target insert length	59	10 kb
Shear speed	Target insert length					
59	10 kb					
	1.2	<p>Note: It is recommended to confirm gDNA is sheared to the appropriate size range (<11 kb) prior to proceeding. If the DNA is under-sheared, a second shear with the same parameters can be repeated. The same Megaruptor consumables can be used if a secondary shear is required.</p>				
	1.3	Transfer sheared DNA into a tube strip or other appropriate tube for the Repair and A tail step. Typical volume loss during shearing is between 5–15 µL.				
		Recommended: evaluate sample quality (concentration and size distribution).				
	1.4	<ul style="list-style-type: none"> Take a 1–2 µL aliquot and measure DNA concentration with a Qubit fluorometer using the 1x dsDNA HS kit. Measure DNA size distribution with a Femto Pulse system. If the DNA size is >11 kb, repeat step 1.2 (see note). 				

Option 2: DNA shearing with a Covaris g-TUBE

✓	Step	Instructions
	1.1	Dilute 1–50 ng of gDNA in a final volume of 65 µL with low TE buffer .
	1.2	Transfer gDNA to the g-TUBE and centrifuge at 2348 x g for 5 minutes to achieve a target mode of 10 kb.
	1.3	Check for any residual sample remaining in the upper chamber of the g-TUBE. If present, re-spin for 1 minute. Repeat spin until the entire gDNA sample has passed through the orifice.
	1.4	Invert and spin the g-TUBE at the same speed selected in step 1.2 until the entire gDNA sample has passed through the g-TUBE orifice.
	1.5	Transfer the recovered sheared DNA to a new 0.2 mL 8-tube strip. Up to 10% volume loss is typical.
		Note: It is acceptable to proceed with up to 60 µL of sample into the Repair and A-tailing step.

SAFE STOPPING POINT - Store at 4°C




Figure 1. Example of 5 ng of human genomic DNA sheared to approximately 10 kb using the Covaris g-TUBE. Distribution of starting genomic DNA is >20 kb.

2. Repair and A-tailing of sheared DNA

This step repairs and prepares the DNA for A-T ligation.

✓ Step Instructions

Add the following components from the SMRTbell prep kit 3.0 to a new microcentrifuge tube. Adjust component volumes for the number of samples being amplified, plus 15% overage. For individual preps, add components directly to the sample from the previous step at the specified volumes and skip steps 2.2 to 2.4.

Repair master mix							
2.1	✓	Tube	Component	Volume			
				Per library	8 libraries*	16 libraries*	24 libraries*
		Purple	Repair buffer	8 μ L	73.6 μ L	147.2 μ L	220.8 μ L
		Blue	End repair mix	2 μ L	18.4 μ L	36.8 μ L	55.2 μ L
		Green	DNA repair mix	1 μ L	9.2 μ L	18.4 μ L	27.6 μ L
			Total volume	11 μL	101 μL	202 μL	304 μL

*15% overage included in master mix calculations

2.2 Pipette-mix the Repair master mix.

2.3 Quick-spin the Repair master mix in a microcentrifuge to collect liquid.

Add 11 μ L of the **Repair master mix** to each sample. Total reaction volume should be around 60 μ L.

2.4 **Note:** It is acceptable to use up to 60 μ L of sample (post-shearing) without negatively impacting enzymatic or cleanup reactions. Recovery sample volume from Megaruptor is typically 50–60 μ L.

2.5 Pipette-mix each sample.

2.6 Quick-spin the sample(s) in a microcentrifuge to collect liquid.

Run the **Repair and A-tailing** thermocycler program with the lid temperature set to >75°C.

	Step	Time	Temperature
2.7	1	30 min	37°C
	2	5 min	65°C
	3	Hold	4°C

2.8 Proceed to the next step of the protocol.

3. Ligation of amplification adapter and cleanup

This step ligates the amplification adapter required for the PCR step.

✓	Step	Instructions																																																
		Add the following components to a new microcentrifuge tube. Adjust component volumes for the number of libraries being prepared, plus 10% overage. For individual preps, add components directly to each sample from the previous step in the order and volume listed below, and skip steps 3.2 to 3.4.																																																
		<table border="1"> <thead> <tr> <th colspan="2">Ligation master mix</th> <th></th> <th></th> <th></th> <th></th> </tr> <tr> <th>✓</th> <th>Tube</th> <th>Component</th> <th>Volume</th> <th></th> <th></th> </tr> <tr> <td></td> <td></td> <td></td> <td></td> <th>Per library</th> <th>8 libraries*</th> <th>16 libraries*</th> <th>24 libraries*</th> </tr> </thead> <tbody> <tr> <td></td> <td></td> <td>Twist universal adapter</td> <td>2 µL</td> <td>17.6 µL</td> <td>35.2 µL</td> <td>52.8 µL</td> </tr> <tr> <td></td> <td>Yellow</td> <td>Ligation mix</td> <td>15 µL</td> <td>132 µL</td> <td>264 µL</td> <td>396 µL</td> </tr> <tr> <td></td> <td>Red</td> <td>Ligation enhancer</td> <td>1 µL</td> <td>8.8 µL</td> <td>17.6 µL</td> <td>26.4 µL</td> </tr> <tr> <td></td> <td></td> <td>Total volume</td> <td>18 µL</td> <td>158.4 µL</td> <td>316.8 µL</td> <td>475.2 µL</td> </tr> </tbody> </table>	Ligation master mix						✓	Tube	Component	Volume							Per library	8 libraries*	16 libraries*	24 libraries*			Twist universal adapter	2 µL	17.6 µL	35.2 µL	52.8 µL		Yellow	Ligation mix	15 µL	132 µL	264 µL	396 µL		Red	Ligation enhancer	1 µL	8.8 µL	17.6 µL	26.4 µL			Total volume	18 µL	158.4 µL	316.8 µL	475.2 µL
Ligation master mix																																																		
✓	Tube	Component	Volume																																															
				Per library	8 libraries*	16 libraries*	24 libraries*																																											
		Twist universal adapter	2 µL	17.6 µL	35.2 µL	52.8 µL																																												
	Yellow	Ligation mix	15 µL	132 µL	264 µL	396 µL																																												
	Red	Ligation enhancer	1 µL	8.8 µL	17.6 µL	26.4 µL																																												
		Total volume	18 µL	158.4 µL	316.8 µL	475.2 µL																																												
3.1																																																		
		*10% overage included in master mix calculation																																																

3.2 Pipette-mix the Ligation master mix.

3.3 Quick-spin the Ligation master mix in a microcentrifuge to collect liquid.

3.4 Add **18 µL** of the **Ligation master mix** containing the **Twist universal adapter** to each sample from the previous step.

The total volume per sample should be around **78 µL**.

3.5 Pipette-mix each sample.

3.6 Quick-spin the sample(s) in a microcentrifuge to collect liquid.

Run the **Adapter ligation** thermocycler program with the lid temperature set to >30°C.

Step	Time	Temperature
1	30 min	20°C
2	Hold	4°C

Cleanup with 1X SMRTbell cleanup beads

3.8 Add **78 µL** of resuspended, **room-temperature SMRTbell cleanup beads** to each sample.

3.9 Pipette-mix the beads until evenly distributed.

3.10 Quick-spin the sample(s) in a microcentrifuge to collect all liquid from the sides of the tubes.

3.11 Incubate at **room temperature** for **10 minutes** to allow DNA to bind beads.

3.12 Place the sample(s) in a magnetic separation rack until beads separate fully from the solution.

3.13 Slowly remove the cleared supernatant without disturbing the beads. Discard the supernatant.

3.14 Slowly dispense **200 µL**, or enough to cover the beads, of **freshly prepared 80% ethanol** into each tube. After **30 seconds**, remove the 80% ethanol and discard.

3.15 Repeat the previous step.

Remove residual 80% ethanol:

- Remove the sample(s) from the magnetic separation rack.
- Quick-spin the sample(s) in a microcentrifuge.
- Place the sample(s) back in a magnetic separation rack until beads separate fully from the solution.
- Remove residual 80% ethanol and discard.

3.16

3.17 Remove the sample(s) from the magnetic rack. **Immediately** add **24 µL of Elution buffer** to each tube and resuspend the beads.

3.18 Quick-spin the sample(s) in a microcentrifuge.

3.19 Incubate at **room temperature** for **5 minutes** to elute DNA.

3.20 Place the sample(s) in a magnetic separation rack until beads separate fully from the solution.

3.21 Slowly aspirate the cleared eluate without disturbing the beads. Transfer eluate to a **new tube**. Discard old tube(s) with beads.

3.22 Proceed to the next step of the protocol.

4. Amplification and cleanup

This step amplifies DNA ligated with the amplification adapter from the previous step.

✓	Step	Instructions																																																
		Add the following components to a new microcentrifuge tube. Adjust component volumes for the number of libraries being prepared, plus 10% overage. Add 4 µL of Twist UDI primers to each sample . For individual preps, add components directly to the sample from the previous step at the specified volumes and skip steps 4.2 to 4.4.																																																
		<table border="1"> <thead> <tr> <th colspan="6">Amplification master mix</th> </tr> <tr> <th>✓</th> <th>Component</th> <th colspan="4">Volume</th> </tr> <tr> <th></th> <th></th> <th>Per library</th> <th>8 libraries*</th> <th>16 libraries*</th> <th>24 libraries*</th> </tr> </thead> <tbody> <tr> <td>4.1</td> <td>2x Xtreme buffer</td> <td>50 µL</td> <td>440 µL</td> <td>880 µL</td> <td>1320 µL</td> </tr> <tr> <td></td> <td>dNTP (2 mM each)</td> <td>20 µL</td> <td>176 µL</td> <td>352 µL</td> <td>528 µL</td> </tr> <tr> <td></td> <td>KOD Xtreme Hot Start DNA polymerase</td> <td>2 µL</td> <td>17.6 µL</td> <td>35.2 µL</td> <td>52.8 µL</td> </tr> <tr> <td></td> <td>Twist UDI Primers (plate)</td> <td>4 µL</td> <td>N/A</td> <td>N/A</td> <td>N/A</td> </tr> <tr> <td></td> <td>Total volume</td> <td>76 µL</td> <td>634 µL</td> <td>1267 µL</td> <td>1901 µL</td> </tr> </tbody> </table> <p>*10% overage included in master mix calculation</p>	Amplification master mix						✓	Component	Volume						Per library	8 libraries*	16 libraries*	24 libraries*	4.1	2x Xtreme buffer	50 µL	440 µL	880 µL	1320 µL		dNTP (2 mM each)	20 µL	176 µL	352 µL	528 µL		KOD Xtreme Hot Start DNA polymerase	2 µL	17.6 µL	35.2 µL	52.8 µL		Twist UDI Primers (plate)	4 µL	N/A	N/A	N/A		Total volume	76 µL	634 µL	1267 µL	1901 µL
Amplification master mix																																																		
✓	Component	Volume																																																
		Per library	8 libraries*	16 libraries*	24 libraries*																																													
4.1	2x Xtreme buffer	50 µL	440 µL	880 µL	1320 µL																																													
	dNTP (2 mM each)	20 µL	176 µL	352 µL	528 µL																																													
	KOD Xtreme Hot Start DNA polymerase	2 µL	17.6 µL	35.2 µL	52.8 µL																																													
	Twist UDI Primers (plate)	4 µL	N/A	N/A	N/A																																													
	Total volume	76 µL	634 µL	1267 µL	1901 µL																																													
	4.2	Pipette-mix the Amplification master mix.																																																
	4.3	Quick-spin the Amplification master mix in a microcentrifuge to collect liquid.																																																
	4.4	On ice, add 72 µL of the Amplification master mix to 28 µL of sample and primer for a total volume of 100 µL.																																																
	4.5	Pipette-mix.																																																
	4.6	Quick-spin the sample(s) in a microcentrifuge to collect liquid.																																																

Run the **PCR** thermocycler program with the lid temperature set to 105°C. Do not add samples to thermal cycler until lid has pre-heated.

Step	Time	Temperature	Cycles
1	2 min	94°C	1 cycle
2	10 sec	98°C	
3	30 sec	58.8°C	8–14 cycles
4	10 min	68°C	
5	7 min	68°C	1 cycle
6	Hold	4°C	

4.7

gDNA input	PCR Cycles*
1 ng	14 cycles
5 ng	12 cycles
10 ng	11 cycles
20 ng	10 cycles
50 ng	8 cycles

*PCR cycles may require optimization

Cleanup with 1X SMRTbell cleanup beads

4.8 Add **100 µL** of resuspended, **room-temperature SMRTbell cleanup beads** to each sample.

4.9 Pipette-mix the beads until evenly distributed.

4.10 Quick-spin the sample(s) in a microcentrifuge to collect liquid.

4.11 Incubate at **room temperature** for **10 minutes** to allow DNA to bind beads.

4.12 Place the sample(s) in a magnetic separation rack until beads separate fully from the solution.

4.13 Slowly remove the cleared supernatant without disturbing the beads. Discard the supernatant.

4.14 Slowly dispense **200 µL**, or enough to cover the beads, of **freshly prepared 80% ethanol** into each tube. After **30 seconds**, remove the 80% ethanol and discard.

4.15 Repeat the previous step.

Remove residual 80% ethanol:

- 4.16
- Remove the sample(s) from the magnetic separation rack.
 - Quick-spin the sample(s) in a microcentrifuge.
 - Place the sample(s) back in a magnetic separation rack until beads separate fully from the solution.

- Remove residual 80% ethanol and discard.

- 4.17 Remove the sample(s) from the magnetic rack. **Immediately** add **50 μ L** of **Elution buffer** to each tube and resuspend the beads by pipetting 10 times or until evenly distributed.
- 4.18 Quick-spin the sample(s) in a microcentrifuge to collect liquid.
- 4.19 Incubate at **room temperature** for **5 minutes** to elute DNA.
- 4.20 Place the sample(s) in a magnetic separation rack until beads separate fully from the solution.
- 4.21 Slowly aspirate the cleared eluate without disturbing the beads. Transfer eluate to a **new tube**. Discard old tube(s) with beads.

Recommended: Evaluate sample quality (concentration and size distribution).

- Take a **1 μ L** aliquot from each tube and dilute with **9 μ L** of **Elution buffer or water**.
 - Measure DNA concentration with a Qubit fluorometer using the 1x dsDNA HS kit.
 - Dilute each aliquot to **250 pg/ μ L** in Femto Pulse dilution buffer.
 - Measure DNA size distribution with a Femto Pulse system.
- 4.22

IMPORTANT: The SMRTbell library preparation requires at least 150 ng of amplified DNA for sequencing on a Revio SMRT Cell using SPRQ chemistry. Sequencing on a Vega or Revio (non-SPRQ) SMRT Cell requires at least 600 ng of amplified DNA. Please see the [Appendix](#) for guidance of additional PCR cycles.




Figure 2. Example of human sample amplified by PCR. Size distribution of amplified products is approximately 10 kb and appropriate to proceed to SMRTbell library construction

- 4.23 Proceed to the next step of the protocol if sample quality and mass is acceptable.

SAFE STOPPING POINT - Store at 4°C

5. Repair and A-tailing of amplified DNA

If multiplexing, indexed amplified samples can be pooled at this stage prior to SMRTbell library preparation with SMRTbell prep kit 3.0 to save on costs and time. Alternatively, samples can be pooled at the end of the protocol prior to sequencing.

This step will prepare the amplified DNA from the previous step for SMRTbell adapter ligation.

✓	Step	Instructions																																																	
		Add the following components from the SMRTbell prep kit 3.0 to a microcentrifuge tube. Adjust component volumes for the number of libraries being prepared, plus 15% overage. For individual preps, add components directly to the sample from the previous step at the specified volumes and skip steps 5.2 to 5.4.																																																	
		<table border="1"> <thead> <tr> <th colspan="7">Repair master mix</th> </tr> <tr> <th>✓</th> <th>Tube</th> <th>Component</th> <th colspan="4">Volume</th> </tr> <tr> <td></td> <td></td> <td></td> <th>Per library</th> <th>8 libraries*</th> <th>16 libraries*</th> <th>24 libraries*</th> </tr> </thead> <tbody> <tr> <td>5.1</td> <td>Purple</td> <td>Repair buffer</td> <td>8 µL</td> <td>73.6 µL</td> <td>147.2 µL</td> <td>220.8 µL</td> </tr> <tr> <td></td> <td>Blue</td> <td>End repair mix</td> <td>2 µL</td> <td>18.4 µL</td> <td>36.8 µL</td> <td>55.2 µL</td> </tr> <tr> <td></td> <td>Green</td> <td>DNA repair mix</td> <td>1 µL</td> <td>9.2 µL</td> <td>18.4 µL</td> <td>27.6 µL</td> </tr> <tr> <td></td> <td></td> <td>Total volume</td> <td>11 µL</td> <td>101 µL</td> <td>202 µL</td> <td>304 µL</td> </tr> </tbody> </table>	Repair master mix							✓	Tube	Component	Volume							Per library	8 libraries*	16 libraries*	24 libraries*	5.1	Purple	Repair buffer	8 µL	73.6 µL	147.2 µL	220.8 µL		Blue	End repair mix	2 µL	18.4 µL	36.8 µL	55.2 µL		Green	DNA repair mix	1 µL	9.2 µL	18.4 µL	27.6 µL			Total volume	11 µL	101 µL	202 µL	304 µL
Repair master mix																																																			
✓	Tube	Component	Volume																																																
			Per library	8 libraries*	16 libraries*	24 libraries*																																													
5.1	Purple	Repair buffer	8 µL	73.6 µL	147.2 µL	220.8 µL																																													
	Blue	End repair mix	2 µL	18.4 µL	36.8 µL	55.2 µL																																													
	Green	DNA repair mix	1 µL	9.2 µL	18.4 µL	27.6 µL																																													
		Total volume	11 µL	101 µL	202 µL	304 µL																																													
		*15% overage included in master mix calculations																																																	
	5.2	Pipette-mix the Repair master mix.																																																	
	5.3	Quick-spin the Repair master mix in a microcentrifuge to collect liquid.																																																	
	5.4	Add 11 µL of the Repair master mix to 49 µL of sample for a total volume of 60 µL.																																																	
	5.5	Pipette-mix each sample.																																																	
	5.6	Quick-spin the sample(s) in a microcentrifuge to collect liquid.																																																	
		Run the Repair and A-tailing thermocycler program with the lid temperature set to >75°C.																																																	
	5.7	<table border="1"> <thead> <tr> <th>Step</th> <th>Time</th> <th>Temperature</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>30 min</td> <td>37°C</td> </tr> <tr> <td>2</td> <td>5 min</td> <td>65°C</td> </tr> <tr> <td>3</td> <td>Hold</td> <td>4°C</td> </tr> </tbody> </table>	Step	Time	Temperature	1	30 min	37°C	2	5 min	65°C	3	Hold	4°C																																					
Step	Time	Temperature																																																	
1	30 min	37°C																																																	
2	5 min	65°C																																																	
3	Hold	4°C																																																	
	5.8	Proceed to the next step of the protocol.																																																	

6. SMRTbell adapter ligation and cleanup

This step ligates the SMRTbell adapter to the ends of each DNA fragment. Samples have already been indexed with the Twist UDIs during the amplification step and therefore do not require the use of a SMRTbell adapter index for demultiplexing. However, if desired, SMRTbell adapter indexes can be used from one of the four SMRTbell adapter index plates 96 (A, B, C or D). The latter will require two separate demultiplexing steps when processing the sequencing data.

✓	Step	Instructions
		Proceed to step 6.2 if not using a SMRTbell adapter index from plates 96 (A, B, C, or D).
6.1	Optional index by SMRTbell adapters:	Add 4 μL of the indexed adapter from the SMRTbell adapter index plate 96 (A, B, C, or D) to each respective sample from the previous step and exclude the SMRTbell adapter from the Ligation master mix (next step). One index per SMRTbell adapter index plate well per sample.

To make the Ligation master mix, add the following components from SMRTbell prep kit 3.0 to a microcentrifuge tube (exclude the SMRTbell adapter if using an indexed adapter from the previous step). Adjust component volumes for the number of libraries being prepared, plus 10% overage. For individual preps, add components directly to each sample from the previous step in the order and volume listed below, and skip steps 6.3 to 6.5.

Ligation master mix						
✓	Tube	Component	Volume			
6.2			Per library	8 libraries**	16 libraries**	24 libraries**
		SMRTbell adapter*	4 μL	35.2 μL	70.4 μL	105.6 μL
	Yellow	Ligation mix	15 μL	132 μL	264 μL	396 μL
	Red	Ligation enhancer	1 μL	8.8 μL	17.6 μL	26.4 μL
		Total volume	20 μL	176 μL	352 μL	528 μL

* Exclude the SMRTbell adapter if using the SMRTbell adapter index plate 96 (A, B, C, or D)

** 10% overage included in master mix calculation.

6.3 Pipette-mix the Ligation master mix.

6.4 Quick-spin the Ligation master mix in a microcentrifuge to collect liquid.

- **Standard, non-indexed SMRTbell adapter:** add 20 μL of the Ligation master mix containing the SMRTbell adapter to each sample.
- **SMRTbell adapter index:** add 16 μL of the Ligation master mix to each sample and SMRTbell adapter mixture.

The total volume per sample should be 80 μL .

6.6 Pipette-mix each sample.

6.7 Quick-spin the sample(s) in a microcentrifuge to collect liquid.

Run the **Adapter ligation** thermocycler program.

	Step	Time	Temperature
6.8	1	30 min	20°C
	2	Hold	4°C

Cleanup with 1X SMRTbell cleanup beads

6.9 Add **80 µL** of resuspended, **room-temperature SMRTbell cleanup beads** to each sample.

6.10 Pipette-mix the beads until evenly distributed.

6.11 Quick-spin the sample(s) in a microcentrifuge to collect all liquid from the sides of the tubes.

6.12 Incubate at **room temperature** for **10 minutes** to allow DNA to bind beads.

6.13 Place the sample(s) in a magnetic separation rack until beads separate fully from the solution.

6.14 Slowly remove the cleared supernatant without disturbing the beads. Discard the supernatant.

6.15 Slowly dispense **200 µL**, or enough to cover the beads, of **freshly prepared 80% ethanol** into each tube. After **30 seconds**, remove the 80% ethanol and discard.

6.16 Repeat the previous step.

Remove residual 80% ethanol:

- Remove the sample(s) from the magnetic separation rack.
- Quick-spin the sample(s) in a microcentrifuge.
- Place the sample(s) back in a magnetic separation rack until beads separate fully from the solution.
- Remove residual 80% ethanol and discard.

6.18 Remove the sample(s) from the magnetic rack. Immediately add **40 µL** of **Elution buffer** to each tube and resuspend the beads.

6.19 Quick-spin the sample(s) in a microcentrifuge.

6.20 Incubate **at room temperature** for **5 minutes** to elute DNA.

6.21 Place the sample(s) in a magnetic separation rack until beads separate fully from the solution.

6.22 Slowly aspirate the cleared eluate without disturbing the beads. Transfer eluate to a **new tube**. Discard old tube(s) with beads.

6.23 Proceed to the next step of the protocol.

SAFE STOPPING POINT - Store at 4°C

7. Nuclease treatment of SMRTbell library

This step removes un-ligated DNA fragments and excess adapter from the library.

✓	Step	Instructions																																										
		Add the following components from the SMRTbell prep kit 3.0 to a new microcentrifuge tube. Adjust component volumes for the number of libraries being prepared, plus 10% overage. For individual preps, add components directly to each sample from the previous step and skip steps 7.2 to 7.4.																																										
		<table border="1"> <thead> <tr> <th colspan="7">Nuclease master mix</th> </tr> <tr> <th>✓</th> <th>Tube</th> <th>Component</th> <th colspan="4">Volume</th> </tr> <tr> <td></td> <td></td> <td></td> <th>Per library</th> <th>8 libraries*</th> <th>16 libraries*</th> <th>24 libraries*</th> </tr> </thead> <tbody> <tr> <td>7.1</td> <td>Light purple</td> <td>Nuclease buffer</td> <td>5 µL</td> <td>44 µL</td> <td>88 µL</td> <td>132 µL</td> </tr> <tr> <td></td> <td>Light green</td> <td>Nuclease mix</td> <td>5 µL</td> <td>44 µL</td> <td>88 µL</td> <td>132 µL</td> </tr> <tr> <td></td> <td></td> <td>Total volume</td> <td>10 µL</td> <td>88 µL</td> <td>176 µL</td> <td>264 µL</td> </tr> </tbody> </table> <p>*10% overage included in master mix calculations</p>	Nuclease master mix							✓	Tube	Component	Volume							Per library	8 libraries*	16 libraries*	24 libraries*	7.1	Light purple	Nuclease buffer	5 µL	44 µL	88 µL	132 µL		Light green	Nuclease mix	5 µL	44 µL	88 µL	132 µL			Total volume	10 µL	88 µL	176 µL	264 µL
Nuclease master mix																																												
✓	Tube	Component	Volume																																									
			Per library	8 libraries*	16 libraries*	24 libraries*																																						
7.1	Light purple	Nuclease buffer	5 µL	44 µL	88 µL	132 µL																																						
	Light green	Nuclease mix	5 µL	44 µL	88 µL	132 µL																																						
		Total volume	10 µL	88 µL	176 µL	264 µL																																						

7.2 Pipette-mix Nuclease master mix.

7.3 Quick-spin the Nuclease master mix in a microcentrifuge to collect liquid.

7.4 Add **10 µL** of **Nuclease master mix** to each sample for a total volume of 50 µL.

7.5 Pipette-mix each sample.

7.6 Quick-spin the sample(s) in a microcentrifuge to collect liquid.

Run the **Nuclease treatment** thermocycler program with the lid temperature set to >47°C.

	Step	Time	Temperature
7.7	1	15 min	37°C
	2	Hold	4°C

7.8 Proceed to the next step of the protocol. It is necessary to perform a bead cleanup on the nuclease reaction prior to safely storing the library or libraries.

8. Diluted AMPure PB bead cleanup and size selection

The AMPure PB bead size selection step progressively depletes DNA fragments less than 5 kb. Ensure accurate ratios are maintained when diluting AMPure PB and when adding the dilution to the library. Failure to do this will result in a loss of sample, or ineffective size-selection.

Important: If DNA size range is <5 kb, proceed with a 1x SMRTbell bead cleanup.

✓	Step	Instructions
	8.1	Make a 35% v/v dilution of AMPure PB beads by adding 350 μ L of resuspended AMPure PB beads to 650 μ L of Elution buffer. The 35% dilution can be stored at 4°C for 30 days. Note: The AMPure PB dilution may be scaled as appropriate for larger scale projects. For accurate pipetting, prepare a minimum volume of 1 mL.
	8.2	Add 3.1X v/v (155 μL) of resuspended, room-temperature 35% AMPure PB beads to each sample from the previous step.
	8.3	Pipette-mix the beads until evenly distributed.
	8.4	Incubate at room temperature for 20 minutes to allow DNA to bind beads.
	8.5	Place sample on an appropriate magnet and allow beads separate fully from the solution.
	8.6	Slowly remove the cleared supernatant without disturbing the beads.
	8.7	Slowly dispense 200 μL , or enough to cover the beads, of freshly prepared 80% ethanol into each sample. After 30 seconds , remove the 80% ethanol and discard.
	8.8	Repeat the previous step.
	8.9	Remove residual 80% ethanol: <ul style="list-style-type: none"> Remove the sample from the magnet. Quick spin to collect liquid at the bottom. Place sample back on the magnet and allow beads separate fully from the solution. Remove residual 80% ethanol and discard.
	8.10	Remove samples from the magnet and immediately add 25 μL of Elution buffer to each sample.
	8.11	Pipette-mix the beads until evenly distributed.
	8.12	Incubate at room temperature for 5 minutes to elute DNA of the beads.
	8.13	Place samples on the magnet and allow the beads to separate fully from the solution.
	8.14	Slowly remove the cleared eluate without disturbing the beads. Transfer eluate to a new tube .

- 8.15 Take a **1 μL** aliquot from each tube and dilute with **9 μL** of **Elution buffer or water**. Measure DNA concentration with a Qubit fluorometer using the 1x dsDNA HS kit. Calculate the total mass. The final overall recovery should be 20-40% as measured from post PCR cleanup to completed SMRTbell library. **DNA concentration must be less than 20 ng/ μL to proceed to ABC.**

Recommended: Further dilute each aliquot to **250 pg/ μL** with Femto Pulse dilution buffer. Measure final SMRTbell library size distribution with a Femto Pulse system.

- 8.16 **If required, dilute 25 μL of library to less than 20 ng/ μL if in the 3–10 kb size range. If DNA size is less than 3 kb, dilute to less than 10 ng/ μL .** If library concentration is higher than recommended for ABC, sequencing performance will be compromised.

8.17 Proceed to Section 9 to prepare library for sequencing with Revio +SPRQ or Vega

Or

Proceed to SMRT Link Sample Setup for preparing samples for Revio non-SPRQ chemistry.

SAFE STOPPING POINT - Store at 4°C for up to 2 weeks or at -20°C long term

9. Annealing, binding, and cleanup (ABC)

This step is for preparing the SMRTbell library (25 μ L) for sequencing on Revio with SPRQ or Vega. If samples are pooled prior to ABC or a custom volume is required, see Appendix A2. **The Polymerase kit used will depend on which sequencer or chemistry is being used (see below).**

Kit	PN
Revio SPRQ polymerase kit	103-496-900
Vega polymerase kit	103-426-500

✓ Step	Instructions																				
9.1	<p>Ensure samples are in the correct concentration range (see 8.16). Prepare the appropriate volume of master mix with 10% overage using the per reaction volumes listed below. For individual preps, add components directly to each sample from the previous step (8.16) and skip steps 9.2 to 9.4.</p> <table border="1"> <thead> <tr> <th colspan="4">Annealing mix</th> </tr> <tr> <th>✓ Tube</th> <th>Component</th> <th>Volume</th> <th></th> </tr> </thead> <tbody> <tr> <td>Light blue</td> <td>Annealing buffer</td> <td>12.5 μL</td> <td></td> </tr> <tr> <td>Light green</td> <td>Standard sequencing primer</td> <td>12.5 μL</td> <td></td> </tr> <tr> <td colspan="3">Total volume</td> <td>25 μL</td> </tr> </tbody> </table>	Annealing mix				✓ Tube	Component	Volume		Light blue	Annealing buffer	12.5 μ L		Light green	Standard sequencing primer	12.5 μ L		Total volume			25 μ L
Annealing mix																					
✓ Tube	Component	Volume																			
Light blue	Annealing buffer	12.5 μ L																			
Light green	Standard sequencing primer	12.5 μ L																			
Total volume			25 μ L																		
9.2	Pipette-mix the Annealing mix and quick-spin to collect liquid.																				
9.3	Add 25 μL of the Annealing mix to each library for a total volume of 50 μL .																				
9.4	Pipette-mix each sample and quick spin.																				
9.5	Incubate at room temperature for 15 minutes .																				
9.6	During primer incubation, prepare the polymerase dilution (see below) and store on ice.																				
9.7	<p>To prepare the polymerase, add the following components to a new microcentrifuge tube on ice. Adjust component volumes for the number of samples being prepared, plus 10% overage.</p> <table border="1"> <thead> <tr> <th colspan="4">Polymerase Dilution</th> </tr> <tr> <th>✓ Tube color</th> <th>Component</th> <th>Volume</th> <th></th> </tr> </thead> <tbody> <tr> <td>Yellow</td> <td>Polymerase buffer</td> <td>47 μL</td> <td></td> </tr> <tr> <td>Purple</td> <td>Sequencing polymerase</td> <td>3 μL</td> <td></td> </tr> <tr> <td colspan="3">Total volume</td> <td>50 μL</td> </tr> </tbody> </table>	Polymerase Dilution				✓ Tube color	Component	Volume		Yellow	Polymerase buffer	47 μ L		Purple	Sequencing polymerase	3 μ L		Total volume			50 μ L
Polymerase Dilution																					
✓ Tube color	Component	Volume																			
Yellow	Polymerase buffer	47 μ L																			
Purple	Sequencing polymerase	3 μ L																			
Total volume			50 μ L																		
9.8	Pipette mix the polymerase dilution and quick-spin to collect liquid.																				
9.9	Add 50 μL of polymerase dilution to primer annealed sample for a total volume of 100 μL .																				
9.10	Pipette-mix each sample and quick-spin to collect liquid.																				
9.11	Incubate at room temperature for 15 minutes .																				

9.12 Proceed immediately to the next step of the protocol to remove excess polymerase.

Post-binding cleanup with 1X SMRTbell cleanup beads

9.13 Add **100 µL** of resuspended, room-temperature SMRTbell cleanup beads to each sample

9.14 Pipette-mix the beads until evenly distributed and quick-spin if necessary to collect all liquid from the sides of the tube.

9.15 Incubate at **room temperature** for **10 minutes** to allow DNA to bind beads

9.16 Place sample on an appropriate magnet and allow beads to separate fully from the solution

9.17 Slowly remove the cleared supernatant without disturbing the beads. Discard the supernatant.
DO NOT USE EtOH. Proceed immediately to the elution. It is important not to let the beads dry out.

Remove sample from the magnet and **immediately** add **Loading buffer** to each tube and resuspend the beads by pipette mixing.

	Revio SPRQ Polymerase Kit	Vega Polymerase Kit
Loading buffer	25 µL	50 µL

9.19 Quick-spin the samples to collect any liquid from the sides of the tube.

9.20 Incubate at **room temperature** for **15 minutes** to elute DNA

9.21 Place sample on magnet and allow beads to separate fully from the solution.

9.22 Slowly aspirate the cleared eluate without disturbing the beads. Transfer eluate to a **new tube**. Discard the old tube with beads

Use **1 µL** of sample to measure DNA concentration with a Qubit fluorometer using the 1x dsDNA HS kit.

9.23 **Important:** The **Qubit Flex** instrument is not compatible with measuring polymerase-bound library in Loading Buffer. Concentration readings will not be accurate.

Proceed to the **Loading Calculator** in SMRT Link v13.3 or higher to calculate the final dilution for adding the sample to Sequencing reagent plate.

9.24 The recommended loading concentration:

- Revio is 120–160 pM.
- Vega is 100–140 pM.

PROTOCOL COMPLETE

Important: Polymerase-bound libraries can be stored at 4°C for up to 1 month, or at -20°C for up to 6 months prior to sequencing. Polymerase-bound libraries can withstand up to 4 freeze-thaw cycles. Note that the Loading buffer is light sensitive.

A1. Appendix –additional amplification for samples with low PCR yield

If there is not enough DNA to proceed with library construction, this section describes a workflow for additional DNA amplification by PCR.

✓ Step Instructions

Add the following components from the KOD Xtreme Hot Start DNA polymerase, plus the Twist UDI primers in the order and volume listed below to a new PCR tube. Adjust component volumes for the number of libraries being prepared, plus 10% overage. For individual preps, add components directly to the sample from the previous step at the specified volumes and skip steps (A1.2 to A1.4).

A1.1	PCR mix	
	Component	Volume
	2x Xtreme buffer	50 μ L
	dNTP (2mM ea.)	20 μ L
	KOD Xtreme Hot Start DNA polymerase	2 μ L
	Twist UDI Primers	4 μ L
	Purified, amplified DNA (step 4.21)	24 μ L
	Total volume	100 μL

A1.2 Pipette-mix the PCR mix. Quick-spin the PCR mix in a microcentrifuge to collect liquid.

A1.3 Mix **24 μ L** purified, amplified DNA with **76 μ L** of the **PCR mix**. Total reaction volume should be 100 μ L.

A1.4 Pipette-mix each sample.

A1.5 Quick-spin the sample(s) in a microcentrifuge to collect liquid.

Run the **PCR** thermocycler program. (Lid: 105°C).

Step	Time	Temperature	Cycles
1	2 min	94°C	1 cycle
2	10 sec	98°C	N* cycles (See below)
3	30 sec	58.8°C	
4	10 min	68°C	
A1.6 5	7 min	68°C	
6	Hold	4°C	

Additional # of cycles	Condition
3	If the total mass <250 ng (<5 ng/ μ L)
4	If the total mass <140 ng (<2.8 ng/ μ L)
5	If the total mass <80 ng (<1.6 ng/ μ L)

Cleanup with 1X SMRTbell cleanup beads

- A1.7 Add **100 μ L** of resuspended, room-temperature SMRTbell cleanup beads to each sample.
- A1.8 Pipette-mix the beads until evenly distributed.
- A1.9 Quick-spin the sample(s) in a microcentrifuge to collect all liquid from the sides of the tubes.
- A1.10 Incubate at **room temperature** for **10 minutes** to allow DNA to bind beads.
- A1.11 Place the sample(s) in a magnetic separation rack until beads separate fully from the solution.
- A1.12 Slowly remove the cleared supernatant without disturbing the beads. Discard the supernatant.
- A1.13 Slowly dispense **200 μ L**, or enough to cover the beads, of **freshly prepared 80% ethanol** into each tube. After 30 seconds, remove the 80% ethanol and discard.
- A1.14 Repeat the previous step.
- Remove residual 80% ethanol:
- Remove the sample(s) from the magnetic separation rack.
 - Quick-spin the sample(s) in a microcentrifuge.
 - Place the sample(s) back in a magnetic separation rack until beads separate fully from the solution.
 - Remove residual 80% ethanol and discard.
- A1.15
- A1.16 Remove the sample(s) from the magnetic rack. Immediately add **50 μ L** of **Elution buffer** to each tube and resuspend the beads.
- A1.17 Quick-spin the sample(s) in a microcentrifuge.
- A1.18 Incubate at **room temperature** for **5 minutes** to elute DNA.
- A1.19 Place the sample(s) in a magnetic separation rack until beads separate fully from the solution.
- A1.20 Slowly remove the cleared eluate without disturbing the beads. Transfer eluate to a **new tube**. Discard old tube(s) with beads.
- A1.21 Use 1 μ L of reamplified samples to quantify with Qubit dsDNA HS kit.

The re-amplified DNA can be stored at 4°C or at -20°C for future use.

A2. Annealing, binding, and cleanup (ABC) for custom volumes

This step is for preparing libraries for sequencing on PacBio HiFi sequencers. **Libraries or pools of libraries must be at a concentration of <20 ng/ μ L. For libraries <3 kb in size, dilute to <10 ng/ μ L.** The sequencing polymerase is stable once bound to the HiFi library and can be stored at 4°C for 1 month or at -20°C for at least 6 months. Use the calculations below to determine reagent volumes based on input sample volume:

	SMRTbell library	Annealing buffer	Standard sequencing primer	Polymerase dilution
Volume (μ L)	x	x/2	x/2	x*2
Example	100	50	50	200

See Section 9 Annealing, binding, and cleanup (ABC) for full protocol.

FAQs

1. Does this protocol work with FFPE samples?

- Ampli-Fi has been successfully executed with FFPE tissue and cell samples and achieves a significantly higher HiFi yield than PCR-free sequencing. A few points to consider are:
 - If not already extracted, we recommend using the QIAamp DNA FFPE Tissue Kit for DNA Extraction and following the [Qiagen Supplementary Protocol](#), which uses Deparaffinization Solution (*cat. No. 19093*).
 - FFPE tissue samples generally do not require shearing. We recommend sizing with Agilent Femto Pulse or TapeStation Genomic DNA ScreenTape. If shearing is required, follow the standard shearing guidelines described in the protocol.
 - For FFPE tissue, the HiFi read length ranges between 800 bp–2 kb in size.
 - For DNA larger than 5 kb, follow the protocol for diluted AMPure size selection (Section 8) for the final SMRTbell cleanup.
 - For DNA between 3–5 kb, perform a 1x SMRTbell bead cleanup for the final SMRTbell cleanup. Diluted AMPure size selection is not recommended and will result in sample loss.
 - For DNA smaller than 3 kb, replace all 1x SMRTbell bead cleanups with 1.3x SMRTbell bead cleanups.
 - To increase the size of fragmented FFPE samples (<3 kb), the post-PCR SMRTbell bead cleanup can be adjusted to be more stringent (1x); however, this will be sample-dependent. If attempting a stringent size-selection, we recommend increasing the number of PCR cycles to account for loss.

2. What should I do if my DNA is under-sheared on the Megaruptor?

- If the Femto Pulse shows a shear mode >11 kb, we recommend performing an additional round of shearing using the same shearing hydropore. Note that further shearing may lead to additional volume loss if a different hydropore is used.
- While KOD Xtreme hot start polymerase can amplify DNA fragments >11 kb, amplification efficiency may be reduced when DNA is under-sheared, resulting in bi-modal HiFi read length distributions.

3. Why do I need to optimize PCR cycles?

- The optimal number of PCR cycles depends on the quality and accurate quantification of the gDNA. Low-quality or low-quantity samples may require additional PCR cycles, while high-quality samples typically require fewer cycles.

4. What are the duplication rates with this protocol?

- Duplication rates vary based on DNA mass input, quality, and the number of PCR cycles used. High-quality hg002 DNA gave the following duplication rates:
 - 12% with 14 PCR cycles
 - 3% with 12 cycles
 - 0.8% with 10 cycles
 - 0.3% with 8 cycles.

For lower duplication rates, we recommend using higher gDNA inputs. If using Revio SPRQ chemistry, PCR cycles may be reduced further.

5. Can I use TapeStation for sizing QC?

- TapeStation requires a higher DNA concentration than what is typically used in the Ampli-Fi protocol, so it is typically not compatible with sheared DNA QC. In addition, we have noted that sheared DNA size is inaccurate and runs significantly larger on TapeStation. However, the SMRTbell product can be quality controlled with Genomic DNA ScreenTape. Note that DNA size may register as larger (up to 7 kb) on the TapeStation.

6. What is the process for demultiplexing Twist UDIs and removing PCR duplicate reads?

- Demultiplexing Twist UDI: if libraries are barcoded using Twist UDIs and *not* the SMRTbell index, treat the sample as “not indexed” in run design:

Samples

Sample is indexed YES NO

Bio Sample Name ?
Required

- In SMRT Link, run the `Demultiplex Barcodes` data utility job, select `Amplifi_TwistUDIadapters_noP7P5` from Barcode Set dropdown, and use `Ampli-Fi_Barcoded_Sample_Name_File.csv` to assign sample names. The `Amplifi_TwistUDIadapters_noP7P5` must be imported in SMRT Link via Data Management before it can be used for demultiplexing. The Barcode Set FASTA and sample assignment CSV are available at <https://www.pacb.com/multiplexing/>. Alternatively, the `lima` command line tool can be used to remove Twist universal adapters and demultiplex Twist UDI primer sequences.
- Removing PCR duplicate read: run the `Mark PCR Duplicates` data utility job in SMRT Link or use `pbmarkdups` command line tool (`pbbioconda`) to eliminate PCR duplicate reads. Additional memory may be required to process the Ampli-Fi samples using the Advanced Parameters. Recommended additional task memory is 16,000 MB.

Revision history (description)	Version	Date
Initial release	01	March 2025
Update to master mix volumes to include 16 and 24 libraries	02	March 2025

Research use only. Not for use in diagnostic procedures. © 2025 Pacific Biosciences of California, Inc. ("PacBio"). All rights reserved. Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at pacb.com/license. Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, SBB, Revio, Onso, Apton, Kinnex, PureTarget, Vega, and SPRQ are trademarks of PacBio.