

# Best Practices for Large Insert Libraries and Sequel System Sequencing

*PacBio UGM*  
*June 27, 2017*

**Melissa Laird Smith, Ph.D.**

Assistant Professor,  
Assistant Director of Technology Development  
Icahn School of Medicine at Mount Sinai  
Dept. of Genetics & Genomics Sciences  
Icahn Institute for Genomics & Multiscale Biology



**Mount  
Sinai**

# Our Technology Development team & capabilities

*Mission: To enable investigators to perform translational and interdisciplinary biomedical research using cutting edge genomics tools*

## Multi-Platform DNA Sequencing:

- 3 PacBio RS IIs
- 2 PacBio Sequels
- 8 Illumina HiSeq 2500/4000
- 2 MiSeq
- 8 Ion Protons
- 7 Ion S5XL
- 2 Ion PGM
- 11 Ion Chefs
- 1 10X Genomics Chromium
- 1 BioNano Genomics Irys
- 1 Applied Biosystems 3730xl
- ONT MinION

## Single Cell Technology:

- 1 10X Genomics Chromium
- 1 BLI Beacon
- 1 BioRad/Illumina ddSeq





# Using long read sequencing as a translational research tool




## Technology Toolbox and Expertise:

- 2 Sequel & 3 RSII sequencing systems with >3000 RSII & >250 Sequel SMRTcells in 2016
- 10X Chromium and BLI Beacon for unique single cell genomics
- BioNano Genomics optical mapping and scaffold generation expertise
- Novel data integration, spanning all technologies to inform R&D and novel diagnostics

## Addressing complex structure & targeting comprehensive diagnostics




## Rapid validation of low freq variants in somatic / heterogeneous tissue



## Infectious Disease & Immune Repertoire Profiling and Surveillance


**Our Goal:** To Address Structural Complexity by  
WGS with Comprehensive Long Read  
Sequencing Methods

# Progress in Human SMRT-WGS Across 7 Genomes



- From previous experience with BNG, can lift scaffold N50 to 22-32 Mb
- Current max contig sizes, even on pass #1 FALCON 0.3 with Sequel are ~10 Mb
- Polished max contigs up to ~20 Mb on RSII completed assemblies
- Speed and cost on Sequel greatly reduced for generating more haplotypes faster
- Additional large genomes prepped for Sequel – mammalian, plants & insects

# Overview of Large Insert Library Prep Process (w/Notes)



**Upfront QC:** request at least 20  $\mu\text{g}$ ; clean (1X AMPure), BA, Qubit

**Shear:** Covaris g-tubes (20 kb shear)

**Pre-library repair:** DR, ER + exoVII always!

**Library Prep:** standard protocols

**Size Selection:** Based on quality & quantity of library, select 10-50 kb or 20-50 kb (BluePippin and/or PippinHT)




**Annealing/Binding:** Primer annealing time (30min); Pol binding time (4h)

**Sequel:** *Titrate* for large projects; 10h movies; 2.0 chemistry






# Sample Data Sets: Input, Library Prep and Size Selection

“Just Fine” sample



“Great” sample




# Sample Data Sets : Sequel Stats

“Just Fine (2 cell pilot)”

## Mapping Report


Value	Analysis Metric
83.52%	Mean Concordance (mapped)
794,400	Number of Subreads (mapped)
5,394,918,161	Number of Subread Bases (mapped)
6,791	Subread Length Mean (mapped)
10,718	Subread Length N50 (mapped)
15,840	Subread Length 95% (mapped)
36,431	Subread Length Max (mapped)
642,617	Number of Polymerase Reads (mapped)
8,555	Polymerase Read Length Mean (mapped)
13,008	Polymerase Read N50 (mapped)
24,010	F
77,546	F



“Great (4 cell pilot)”

## Mapping Report

Value	Analysis Metric
84.04%	Mean Concordance (mapped)
2,172,791	Number of Subreads (mapped)
18,590,421,740	Number of Subread Bases (mapped)
8,556	Subread Length Mean (mapped)
13,904	Subread Length N50 (mapped)
21,620	Subread Length 95% (mapped)
50,116	Subread Length Max (mapped)
1,758,655	Number of Polymerase Reads (mapped)
10,732	Polymerase Read Length Mean (mapped)
17,498	Polymerase Read N50 (mapped)
32,120	F
75,133	F



**(7) Cau-haplotype (SQL-10 hr)**

~2 Mb\* de novo PB-FLCN

~50 SMRTcells


~ weeks

~ 17.5kb subRL N50



# Looking Forward

- Awaiting SMRTLink 5.0 and IPS upgrades soon
- Currently:
  - 5 whole human genomes in progress/in queue
  - Multiple additional large genomes, including mammalian, insect and plants
- Expanding to MegaRuptor ( $\geq 30$  kb)-derived libraries



# Many Thanks!

## Tech Dev Team

**Bobby Sebra**

**Gintaras Deikus**

**Irina Oussenko**

**Hardik Shah**

**Mondale Wang**

**Kimaada Allette**

**Maya Strahl**

**Ethan Ellis**

**Matthew Emery**

**Diane Castillo**

**Leah Newman**

**James Powell**

**Alesia Antoine**

**Wissam Hamou**

**Hanane Arib**

**David Sachs**



**Open to all the Q&A!**

*melissa.smith@mssm.edu*